**Test Definition: APCRR** 



BORATORIES Activated Protein C Resistance V, with Reflex to Factor V Leiden, Blood and Plasma

# **Reporting Title:** APCRV, w/Reflex, P Performing Location: Rochester

# Specimen Requirements:

Blood and plasma are required.

Patient Preparation: Fasting preferred

Specimen Type: Whole blood Container/Tube: Preferred: Lavender top (EDTA) Acceptable: Yellow top (ACD) or light-blue top (3.2% sodium citrate) Specimen Volume: 3 mL Collection Instructions: 1. Invert several times to mix blood. 2. Send whole blood specimen in original tube. Do not aliquot.

Specimen Type: Platelet-poor plasma Collection Container/Tube: Light-blue top (3.2% sodium citrate) Submission Container/Tube: Plastic vial, polypropylene preferred Specimen Volume: 1 mL **Collection Instructions:** 

1. For complete instructions, see Coagulation Guidelines for Specimen Handling and Processing.

2. Within 4 hours of collection, centrifuge, transfer all plasma into a plastic vial, and centrifuge plasma again.

3. Aliquot plasma into separate plastic vial leaving 0.25 mL in the bottom of centrifuged vial.

4. Freeze plasma aliquot immediately at -20 degrees C, or ideally, at -40 degrees C or below. Additional Information:

1. Double-centrifuged specimen is critical for accurate results as platelet contamination may cause spurious results.

2. Each coagulation assay requested should have its own vial.

#### Specimen Minimum Volume:

Plasma: 0.5 mL; Whole blood: 3 mL

# Forms:

1. New York Clients-Informed consent is required. Document on the request form or electronic order that a copy is on file. The following documents are available:

-Informed Consent for Genetic Testing (T576)

-Informed Consent for Genetic Testing-Spanish (T826)

2. Coagulation Patient Information (T675)

3. If not ordering electronically, complete, print, and send a Coagulation Test Request (T753) with the specimen.

 
 MAYOCLINIC
 Test Definition: APCRR

 LABORATORIES
 Activated Protein C Resistance V, with Reflex to Factor V Leiden,
**Blood and Plasma** 

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	14 days	
Whole blood	Ambient (preferred)	14 days	
	Frozen	14 days	
	Refrigerated	14 days	

# **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
APCR	APCRV Ratio	Numeric		13590-5
	Also used by tests: APCRV			
INT55	Interpretation	Alphanumeric		48591-2
	Also used by tests: APCRV			
SC018	Whole Blood	Alphanumeric		No LOINC Needed

LOINC and CPT codes are provided by the performing laboratory.

### **Supplemental Report:**

No

#### **Components:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
APCRV	Activated Protein Resistance V, P			Yes	Yes
SC018	Whole Blood			Yes	No

# **CPT Code Information:**

85307



 
 MAYOCLINIC
 Test Definition: APCKK

 LABORATORIES
 Activated Protein C Resistance V, with Reflex to Factor V Leiden, Blood and Plasma
**Blood and Plasma** 

# **Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
F5DNA	Factor V Leiden (R506Q) Mutation, B			No	Yes
F5DNI	APCRV/F5DNA Summary Interpretation			No	No

# **Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
F5DNA	21838	Factor V Leiden (R506Q) Mutation, B	Alphanumeric		21668-9
F5DNA	21839	F5DNA Interpretation	Alphanumeric		69049-5
F5DNA	21841	F5DNA Reviewed By	Alphanumeric		18771-6
F5DNI	38238	APCRV/F5DNA Summary Interpretation	Alphanumeric		48591-2
F5DNI	44181	Interpretation	Alphanumeric		69049-5

#### **Reference Values:**

ACTIVATED PROTEIN C RESISTANCE V RATIO

> or =2.3

Pediatric reference range has neither been established nor is available in scientific literature. The adult reference range likely would be applicable to children older than 6 months.