

**Reporting Title:** Oxalate, Dialysate Fluid**Performing Location:** Rochester**Specimen Requirements:**

Specimen Type: Dialysate fluid

Patient Preparation: Patient should avoid taking vitamin C supplements for 24 hours prior to dialysis

Supplies: Urine Tubes, 10 mL (T068)

Container/Tube: Plastic, 10-mL urine tube or a clean, plastic aliquot container with no metal cap or glued insert

Specimen Volume: 5 mL

Collection Instructions: Adjust the pH of the specimen to 2.5 to 3.0 with 6M Hydrochloric Acid

Additional Information: Nonacidified frozen hemodialysate delivered to the laboratory within 3 days from collection will be accepted and the following comment will be added to the result: In nonacidified hemodialysate stored frozen, oxalate values may increase spontaneously.

**Specimen Minimum Volume:**

2 mL

**Forms:**

If not ordering electronically, complete, print, and send a Renal Diagnostics Test Request (T830) with the specimen.

Specimen Type	Temperature	Time	Special Container
Dialysate Fluid	Frozen	14 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
DOXA1	Oxalate, Dialysate Fluid	Numeric	mcmol/L	47715-8

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

83945

---

**Reference Values:**

Not applicable