

Triglycerides, Body Fluid

Reporting Title: Triglycerides, BF **Performing Location:** Rochester

Ordering Guidance:

For help distinguishing between chylous and nonchylous effusions, order BFLA1 / Lipid Analysis, Body Fluid. The body fluid will be tested for cholesterol and triglyceride concentrations and undergo lipoprotein electrophoresis.

Necessary Information:

- 1. Date and time of collection is required.
- 2. Specimen source is required.

Specimen Requirements:

Preferred Source: -Peritoneal fluid (peritoneal, abdominal, ascites, paracentesis) -Pleural fluid (pleural, chest, thoracentesis) -Drain fluid (drainage, JP drain) -Pericardial Fluid Acceptable Source: Write in source name with source location (if appropriate) Collection Container/Tube: Sterile container Submission Container/Tube: Plastic vial Specimen Volume: 1 mL Collection Instructions: 1. Centrifuge to remove any cellular material and transfer into a plastic vial. 2. Indicate the specimen source and source location on label.

Specimen Minimum Volume:

0.5 mL

Specimen Type	Temperature	Time	Special Container
Body Fluid	Frozen (preferred)	30 days	
	Refrigerated	7 days	
	Ambient	24 hours	



Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
TGLBF	FLD25	Fluid Type:	Plain Text	Yes

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
FLD25	Fluid Type:	Alphanumeric		14725-6
	Also used by tests: CHLBF			
BFTGL	Triglycerides (BF)	Numeric	mg/dL	12228-3

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

84478

Reference Values:

An interpretive report will be provided