

Test Definition: 11INE

Factor XI Inhibitor Evaluation, Plasma

Reporting Title: Factor XI Inhib Profile, P

Performing Location: Rochester

Ordering Guidance:

This test is for factor XI inhibitors only. If the presence or type of inhibitor is unknown, first order APROL / Prolonged Clot Time Profile, Plasma, except for patients with known hemophilia A or B. When screening studies are needed for patients with known hemophilia A or B, order 8INHE / Factor IX Inhibitor Evaluation, Plasma; or 9 INHE / Factor IX Inhibitor Evaluation, Plasma; respectively.

Shipping Instructions:

Send all vials in the same shipping container.

Necessary Information:

If priority specimen, mark request form, give reason, and request a call-back.

Specimen Requirements:

Specimen Type: Platelet-poor plasma

Patient Preparation:

- 1. Patient should not be receiving Coumadin (warfarin), heparin, direct thrombin inhibitors (argatroban, dabigatran), or direct factor Xa inhibitors (apixaban, rivaroxaban, and edoxaban).
- 2. Fasting preferred.

Collection Container/Tube: Light-blue top (3.2% sodium citrate)

Submission Container/Tube: 3 Plastic vials

Specimen Volume: 3 mL in 3 plastic vials, each containing 1 mL

Collection Instructions:

- 1. Specimen must be collected prior to factor replacement therapy.
- 2. For complete instructions, see Coagulation Guidelines for Specimen Handling and Processing.
- 3. Centrifuge, transfer all plasma into a plastic vial, and centrifuge plasma again.
- 4. Aliquot plasma (1-2 mL per aliquot) into 3 separate plastic vials, leaving 0.25 mL in the bottom of centrifuged vial.
- 5. Freeze plasma immediately at -20 degrees C or, ideally, at -40 degrees C or below.

Additional Information:

- 1. Double-centrifuged specimen is critical for accurate results as platelet contamination may cause spurious results.
- 2. Each coagulation assay requested should have its own vial.

Specimen Minimum Volume:

2 Plastic vials, each containing 1 mL

Forms:

If not ordering electronically, complete, print, and send a Coagulation Test Request (T753) with the specimen.

Test Definition: 11INE

Factor XI Inhibitor Evaluation, Plasma

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	14 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
11INT	FXI Inhib Profile Tech Interp	Alphanumeric		69049-5
F_11	Coag Factor XI Assay, P	Numeric	%	3226-8
	Also used by tests: F_11			

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Components:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
11INT	FXI Inhib Profile Tech Interp			Yes	No
F_11	Coag Factor XI Assay, P			Yes	Yes

CPT Code Information:

85390

85270

85335 (if appropriate)

85335 (if appropriate)

85390 (if appropriate)

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
11AIH	FXI Inhib Profile Prof Interp			No	No

Test Definition: 11INE

Factor XI Inhibitor Evaluation, Plasma

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
11_IS	Factor XI Inhib Scrn			No	No
GBETH	General Factor Bethesda Units, P			No	No

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
11AIH	607492	Reviewed by	Alphanumeric		18771-6
11AIH	607448	FXI Inhib Profile Prof Interp	Alphanumeric		69049-5
11_IS	7804	Factor XI Inhib Scrn	Alphanumeric		80603-4
GBETH	607434	General Factor Bethesda Units, P	Numeric	BU	13591-3

Reference Values:

FACTOR XI ACTIVITY ASSAY

Adults: 55-150%

Normal, full-term newborn infants or healthy premature infants may have decreased levels (> or =10%) that may not reach adult levels for 180 days or more postnatal.*

*See Pediatric Hemostasis References section in Coagulation Guidelines for Specimen Handling and Processing.

FACTOR XI INHIBITOR SCREEN:

Negative

GENERAL FACTOR BETHESDA UNITS:

< or =0.5 Bethesda Units