

Test Definition: SPAS

Arsenic Speciation, Random, Urine

Reporting Title: Arsenic Speciation, Random, U

Performing Location: Rochester

Specimen Requirements:

Patient Preparation:

1. Patient should not eat seafood for a 48-hour period prior to start of collection.

2. High concentrations of gadolinium and iodine are known to interfere with most metal tests. If either gadolinium- or iodine-containing contrast media has been administered, a specimen should not be collected for at least 96 hours. Supplies:

-Urine Tubes, 10 mL (T068)

-Sarstedt Aliquot Tube 5 mL (T914)

Collection Container/Tube: Clean, plastic urine collection container

Submission Container/Tube: Plastic vial or a clean, plastic aliquot container with no metal cap or glued insert

Specimen Volume: 5 mL Collection Instructions:

1. Collect a random urine specimen.

2. See Metals Analysis Specimen Collection and Transport for complete instructions.

Specimen Minimum Volume:

3 mL

Specimen Type	Temperature	Time	Special Container	
Urine	Refrigerated (preferred)	28 days		
	Frozen	28 days		
	Ambient	72 hours		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
607692	Inorganic Arsenic (Toxic)	Alphanumeric	mcg/L	12481-8
607693	Methylated Arsenic (Toxic)	Alphanumeric	mcg/L	53779-5
607694	Organic Arsenic (Non-Toxic)	Alphanumeric	mcg/L	53778-7
609405	Toxic Arsenic Concentration	Alphanumeric	mcg/L	96250-6
607691	Total Arsenic Concentration	Numeric	mcg/L	5586-3

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:



Test Definition: SPAS

Arsenic Speciation, Random, Urine

No

CPT Code Information:

82175

Reference Values:

TOXIC ARSENIC <35 mcg/L
Reference values apply to all ages.

Arsenic Speciation Interpretive Information:

The toxic arsenic concentration represents the sum of the inorganic and methylated arsenic species. The reference value for toxic arsenic is <35 mcg/L. This value is based on the ACGIH Biological Exposure Index (BEI), which does not include the nontoxic organic arsenic.