

Test Definition: 2B6Q

Cytochrome P450 2B6 Genotype, Varies

Reporting Title: CYP2B6 Genotype, V

Performing Location: Rochester

Ordering Guidance:

Testing is available as the single gene assay (this test) and as a part of a psychotropic pharmacogenomics panel.

If genotype testing for psychotropic medications is requested, order PSYQP / Psychotropic Pharmacogenomics Gene Panel, Varies.

Specimen Requirements:

Submit only 1 of the following specimens:

Specimen Type: Whole blood

Container/Tube: Lavender top (EDTA)

Specimen Volume: 3 mL Collection Instructions:

1. Invert several times to mix blood.

2. Send whole blood specimen in original tube. Do not aliquot.

Specimen Stability Information: Ambient (preferred) 9 days/Refrigerated 30 days

Specimen Type: Saliva

Patient Preparation: Patient should not eat, drink, smoke, or chew gum 30 minutes prior to collection.

Supplies: Saliva Swab Collection Kit (T786)

Specimen Volume: 1 Swab

Collection Instructions: Collect and send specimen per kit instructions.

Specimen Stability Information: Ambient 30 days

Specimen Type: Extracted DNA Container/Tube: 2 mL screw top tube Specimen Volume: 100 mcL (microliters)

Collection Instructions:

- 1. The preferred volume is 100 mcL at a concentration of 50 ng/mcL.
- Provide concentration of DNA and volume on tube.

Specimen Stability Information: Frozen (preferred) 1 year/Ambient/Refrigerated

Specimen Minimum Volume:

Blood: 1 mL

Saliva, extracted DNA: see Specimen Required

Forms:

- 1. New York Clients-Informed consent is required. Document on the request form or electronic order that a copy is on file. The following documents are available:
- -Informed Consent for Genetic Testing (T576)
- -Informed Consent for Genetic Testing-Spanish (T826)

Test Definition: 2B6Q

Cytochrome P450 2B6 Genotype, Varies

2. If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

| Specimen Type | Temperature | Time | Special Container | |
|---------------|-------------|------|-------------------|--|
| Varies | Varies | | | |

Result Codes:

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|------------------------|--------------|------|------------|
| 610082 | CYP2B6 Genotype | Alphanumeric | | 72882-4 |
| 610083 | CYP2B6 Phenotype | Alphanumeric | | 79720-9 |
| 610566 | CYP2B6 Activity Score | Alphanumeric | | In Process |
| 610084 | Interpretation | Alphanumeric | | 69047-9 |
| 610085 | Additional Information | Alphanumeric | | 48767-8 |
| 610086 | Method | Alphanumeric | | 85069-3 |
| 610087 | Disclaimer | Alphanumeric | | 62364-5 |
| 610088 | Reviewed by | Alphanumeric | | 18771-6 |

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

81479

Reference Values:

An interpretive report will be provided.