

# **Test Definition: DPYDQ**

Dihydropyrimidine Dehydrogenase Genotype, Varies

Reporting Title: DPYD Genotype, V

Performing Location: Rochester

### Ordering Guidance:

This test does not detect or report variants other than the \*2A, \*7, \*8, \*10, \*13, rs67376798, rs75017182, and rs115232898 alleles. Sequencing of the full gene is available for detection of additional variants as well as the alleles listed: order DPYDZ / Dihydropyrimidine Dehydrogenase, DPYD Full Gene Sequencing, Varies.

### **Specimen Requirements:**

Multiple genotype tests can be performed on a single specimen after a single extraction. See Multiple Genotype Test List for a list of tests that can be ordered together.

Submit only 1 of the following specimens:

Specimen Type: Whole blood

Container/Tube: Lavender top (EDTA)

Specimen Volume: 3 mL Collection Instructions:

1. Invert several times to mix blood.

2. Send whole blood specimen in original tube. Do not aliquot.

Specimen Stability Information: Ambient (preferred) 9 days/Refrigerated 30 days

Specimen Type: Saliva

Patient Preparation: Patient should not eat, drink, smoke, or chew gum 30 minutes prior to collection.

Supplies: Saliva Swab Collection Kit (T786)

Specimen Volume: 1 Swab

Collection Instructions: Collect and send specimen per kit instructions.

Specimen Stability Information: Ambient 30 days

Specimen Type: Extracted DNA Container/Tube: 2-mL Screw top tube Specimen Volume: 100 mcL (microliters)

Collection Instructions:

- 1. The preferred volume is 100 mcL at a concentration of 50 ng/mcL.
- 2. Provide concentration of DNA and volume on tube.

Specimen Stability Information: Frozen (preferred)/Ambient/Refrigerated

### Specimen Minimum Volume:

Blood: 0.4 mL; Saliva, extracted DNA: see Specimen Required

### Forms:

- 1. New York Clients-Informed consent is required. Document on the request form or electronic order that a copy is on file. The following documents are available:
- -Informed Consent for Genetic Testing (T576)



# **Test Definition: DPYDQ**

Dihydropyrimidine Dehydrogenase Genotype, Varies

- -Informed Consent for Genetic Testing (Spanish) (T826)
- 2. If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
610138	DPYD Phenotype	Alphanumeric		79719-1
610139	DPYD Activity Score	Numeric		In Process
613999	DPYD Genotype	Alphanumeric		45284-7
610140	Interpretation	Alphanumeric		69047-9
610141	Additional Information	Alphanumeric		48767-8
610142	Method	Alphanumeric		85069-3
610143	Disclaimer	Alphanumeric		62364-5
610144	Reviewed by	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

## **CPT Code Information:**

81232

#### **Reference Values:**

DPYD Phenotype: Normal metabolizer

DPYD Activity Score: 2.00

DPYD Genotype: No variants were detected in the DPYD gene.

An interpretive report will be provided.