

## **Test Definition: TPNUQ**

Thiopurine Methyltransferase (TPMT) and Nudix Hydrolase (NUDT15) Genotyping, Varies

Reporting Title: TPMT and NUDT15 Genotype, V

Performing Location: Rochester

#### Ordering Guidance:

For thiopurine methyltransferase (TPMT) enzyme activity testing, order TPMT3 / Thiopurine Methyltransferase Activity Profile, Erythrocytes; however, this test should also be ordered because TPMT enzyme activity testing cannot detect variants in NUDT15, which also impact thiopurine metabolism.

#### **Specimen Requirements:**

Multiple genotype tests can be performed on a single specimen after a single extraction. See Multiple Genotype Test List for a list of tests that can be ordered together.

Submit only 1 of the following specimens:

Specimen Type: Whole blood

Container/Tube: Lavender top (EDTA)

Specimen Volume: 3 mL Collection Instructions:

1. Invert several times to mix blood.

2. Send whole blood specimen in original tube. Do not aliquot.

Specimen Stability Information: Ambient (preferred) 9 days/Refrigerated 30 days

Specimen Type: Saliva

Patient Preparation: Patient should not eat, drink, smoke, or chew gum 30 minutes prior to collection.

Supplies: Saliva Swab Collection Kit (T786)

Specimen Volume: 1 Swab

Collection Instructions: Collect and send specimen per kit instructions.

Specimen Stability Information: Ambient 30 days

Specimen Type: Extracted DNA Container/Tube: 2-mL screw top tube

Specimen Volume: 100 mcL Collection Instructions:

1. The preferred volume is 100 mcL at a concentration of 50 ng/mcL.

Provide concentration of DNA and volume on tube.

Specimen Stability Information: Frozen (preferred)/Ambient/Refrigerated

#### **Specimen Minimum Volume:**

Blood: 0.4 mL

Saliva, extracted DNA: see Specimen Required

#### Forms:

1. New York Clients-Informed consent is required. Document on the request form or electronic order that a copy is on file. The following documents are available:



## **Test Definition: TPNUQ**

Thiopurine Methyltransferase (TPMT) and Nudix Hydrolase (NUDT15) Genotyping, Varies

- -Informed Consent for Genetic Testing (T576)
- -Informed Consent for Genetic Testing-Spanish (T826)
- 2. If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:
- -Neurology Specialty Testing Client Test Request (T732)
- -Gastroenterology and Hepatology Test Request (T728)
- -Therapeutics Test Request (T831)

| Specimen Type | Temperature | Time | Special Container |  |
|---------------|-------------|------|-------------------|--|
| Varies        | Varies      |      |                   |  |

#### **Result Codes:**

| Result ID | Reporting Name         | Туре         | Unit | LOINC®  |
|-----------|------------------------|--------------|------|---------|
| 610159    | TPMT Genotype          | Alphanumeric |      | 41048-0 |
| 610160    | TPMT Phenotype         | Alphanumeric |      | 79713-4 |
| 610161    | NUDT15 Genotype        | Alphanumeric |      | 93194-9 |
| 610162    | NUDT15 Phenotype       | Alphanumeric |      | 93195-6 |
| 610163    | Interpretation         | Alphanumeric |      | 69047-9 |
| 610164    | Additional Information | Alphanumeric |      | 48767-8 |
| 610165    | Method                 | Alphanumeric |      | 85069-3 |
| 610166    | Disclaimer             | Alphanumeric |      | 62364-5 |
| 610167    | Reviewed by            | Alphanumeric |      | 18771-6 |

LOINC and CPT codes are provided by the performing laboratory.

### **Supplemental Report:**

Νo

#### **CPT Code Information:**

0034U



# **Test Definition: TPNUQ**

Thiopurine Methyltransferase (TPMT) and Nudix Hydrolase (NUDT15) Genotyping, Varies

#### **Reference Values:**

An interpretive report will be provided.