

# **Test Definition: UGTFZ**

UDP-Glucuronosyltransferase 1A1 (UGT1A1), Full Gene Sequencing, Varies

Reporting Title: UGT1A1 Full Gene Sequencing, V

Performing Location: Rochester

### **Ordering Guidance:**

If analysis of only the UGT1A1 promoter TA repeat region (\*28, \*36, \*37 alleles) is desired, see U1A1Q / UDP-Glucuronosyltransferase 1A1 TA Repeat Genotype, UGT1A1, Varies.

## **Shipping Instructions:**

If submitting microtube, place inside a larger tube or vial for transport.

# **Specimen Requirements:**

Patient Preparation: A previous liver transplant, bone marrow transplant from an allogenic donor, or a recent (ie, <6 weeks from time of sample collection) heterologous blood transfusion will interfere with testing. Call 800-533-1710 for instructions for testing patients who have received a bone marrow transplant.

Submit only 1 of the following specimens:

Specimen Type: Whole blood

Container/Tube:

Adults: Lavender top (EDTA) Pediatrics: Purple microtube

Specimen Volume: Adults: 3 mL Pediatrics: 1 mL Collection Instructions:

1. Invert several times to mix blood.

2. Send whole blood specimen in original tube. Do not aliquot.

Specimen Stability Information: Ambient (preferred) 9 days/Refrigerated 30 days

Specimen Type: Saliva

Patient Preparation: Patient should not eat, drink, smoke, or chew gum 30 minutes prior to collection.

Supplies: Saliva Swab Collection Kit (T786) Container/Tube: Saliva Swab Collection Kit

Specimen Volume: One swab

Collection Instructions: Collect and send specimen per kit instructions.

Specimen Stability Information: Ambient 30 days

#### **Specimen Minimum Volume:**

Whole blood: 0.45 mL

Saliva: See Specimen Required

#### Forms:



# **Test Definition: UGTFZ**

UDP-Glucuronosyltransferase 1A1 (UGT1A1), Full Gene Sequencing, Varies

- 1. New York Clients-Informed consent is required. Document on the request form or electronic order that a copy is on file. The following documents are available:
- -Informed Consent for Genetic Testing (T576)
- -Informed Consent for Genetic Testing-Spanish (T826)
- 2. UGT1A1 Gene Testing Patient Information (T664) is requested but not required.
- 3. If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:
- -Renal Diagnostics Test Request (T830)
- -Therapeutics Test Request (T831)

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
618686	Result Summary	Alphanumeric		50397-9
618687	TA Repeat Result	Alphanumeric		95143-4
618688	Full Gene Sequence Result	Alphanumeric		82939-0
618691	Interpretation	Alphanumeric		69047-9
618692	Method	Alphanumeric		85069-3
618693	Disclaimer	Alphanumeric		62364-5
618694	Additional Information	Alphanumeric		48767-8
618695	Reviewed By	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

### **CPT Code Information:**

81404



# **Test Definition: UGTFZ**

UDP-Glucuronosyltransferase 1A1 (UGT1A1), Full Gene Sequencing, Varies

### **Reference Values:**

An interpretive report will be provided.