

Test Definition: KRASP

KRAS Somatic Mutation Analysis, Tumor

Reporting Title: KRAS Somatic Mutation Analysis **Performing Location:** Rochester

Necessary Information:

Pathology report (final or preliminary) must accompany specimen in order for testing to be performed. At minimum, it should contain the following information:

- 1. Patient name
- 2. Block number-must be on all blocks, slides and paperwork (can be handwritten on the paperwork)
- 3. Tissue collection date
- 4. Source of the tissue

Specimen Requirements:

Preferred: Specimen Type: Tissue Container/Tube: Tissue block Collection Instructions: Submit a formalin-fixed, paraffin-embedded tissue block.

Acceptable: Specimen Type: Tissue Container/Tube: Slides Specimen Volume: 1 stained and 5 unstained Collection Instructions: Submit 1 slide stained with hematoxylin and eosin and 5 unstained, nonbaked slides with 5-micron thick sections of the tumor tissue.

Specimen Type: Cytology slide (direct smears or ThinPrep) Slides: 1 to 3 slides Collection Instructions: Submit 1 to 3 slides stained and coverslipped with a preferred total of 5000 nucleated cells or a minimum of at least 3000 nucleated cells.

Note: Glass coverslips are preferred; plastic coverslips are acceptable but will result in longer turnaround times. Additional Information: Cytology slides will not be returned.

Specimen Minimum Volume:

Formalin-fixed, paraffin-embedded (FFPE) tissue block (preferred) or 1 slide stained with hematoxylin and eosin and 5 unstained, nonbaked slides (5-microns thick sections) of the tumor tissue.

Forms:

- 1. Molecular Genetics: Inherited Cancer Syndromes Patient Information (T519) in Special Instructions.
- 2. If not ordering electronically, complete, print, and send a Oncology Test Request (T729) with the specimen.

Specimen Type T	Temperature	Time	Special Container
-----------------	-------------	------	-------------------

LABORATORIES

Test Definition: KRASP

KRAS Somatic Mutation Analysis, Tumor

Varies	Ambient (preferred)	
	Refrigerated	

Result Codes:

Result ID	Reporting Name	Type Unit		LOINC®
610680	Result Summary	Alphanumeric		50397-9
610681	Result	Alphanumeric		82939-0
610682	Interpretation	Alphanumeric		69047-9
610683	Specimen	Alphanumeric		31208-2
610684	Source	Alphanumeric		31208-2
610685	Tissue ID	Alphanumeric		80398-1
610686	Released By	Alphanumeric		18771-6
610687	Method	Alphanumeric		85069-3
610688	Disclaimer	Alphanumeric		62364-5
614165	Additional Information	Alphanumeric		48767-8

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

81275-KRAS (v-Ki-ras2 Kirsten rat sarcoma viral oncogene) (eg, carcinoma) gene analysis, variants in codons 12 and 13 81276-KRAS additional variant(s) 82221 Migradiagastian, manual

88381-Microdissection, manual

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
SLIRV	Slide Review in MG			Yes	No (Bill Only)



KRAS Somatic Mutation Analysis, Tumor

Reference Values:

An interpretive report will be provided.