

# **Test Definition: LPLFX**

ABORATORIES Lymphoplasmacytic Lymphoma/Waldenstrom Macroglobulinemia,
MYD88 L265P with Reflex to CXCR4, Varies

Reporting Title: Reflex Testing of MYD88 and CXCR4

Performing Location: Rochester

# **Shipping Instructions:**

Whole blood or bone marrow specimens must arrive within 10 days of collection.

# **Necessary Information:**

The following information is required:

- 1. Pertinent clinical history
- 2. Clinical or morphologic suspicion
- 3. Date and time of collection
- 4. Specimen source

# Specimen Requirements:

Submit only 1 of the following specimens:

Preferred:

Specimen Type: Bone marrow

Container/Tube:

Preferred: Lavender top (EDTA) Acceptable: Yellow top (ACD) Specimen Volume: 2 mL Collection Instructions:

- 1. Invert several times to mix bone marrow
- 2. Send bone marrow specimen in original tube. Do not aliquot.
- 3. Label specimen as bone marrow

Specimen Stability Information: Ambient (preferred)/Refrigerated

Specimen Type: Paraffin-embedded tissue

Container/Tube: Paraffin block Specimen Stability: Ambient

Specimen Type: Paraffin-embedded bone marrow aspirate clot

Container/Tube: Paraffin block Specimen Stability: Ambient

Specimen Type: Tissue Slides: Unstained slides

Specimen Volume: 10 to 20 slides

Additional Information: Tissue must demonstrate involvement by a hematologic neoplasm (eg, acute myelocytic

leukemia), not solid tumors.

Specimen Stability Information: Ambient

Acceptable:

Specimen Type: Whole blood



# **Test Definition: LPLFX**

ABORATORIES Lymphoplasmacytic Lymphoma/Waldenstrom Macroglobulinemia,
MYD88 L265P with Reflex to CXCR4, Varies

Container/Tube:

Preferred: Lavender top (EDTA) Acceptable: Yellow top (ACD) Specimen Volume: 3 mL Collection Instructions:

1. Invert several times to mix blood

2. Send whole blood specimen in original tube. Do not aliquot.

3. Label specimen as blood

Specimen Stability Information: Ambient (preferred)/Refrigerated

Specimen Type: Extracted DNA

Container/Tube: 1.5- to 2-mL tube with indication of volume and concentration of the DNA

Specimen Volume: Entire specimen

Collection Instructions: Label specimen as extracted DNA and list the specimen source. Include indication of volume and

concentration of the DNA.

Specimen Stability Information: Frozen (preferred)/Refrigerated/Ambient

# Specimen Minimum Volume:

Whole blood, Bone marrow: 1 mL

Extracted DNA: 50 mcL with a concentration of at least 20 nanograms per mcL

Other specimen types: See Specimen Required

## Forms:

1. Hematopathology Patient Information (T676)

2. If not ordering electronically, complete, print, and send a Hematopathology/Cytogenetics Test Request (T726) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies	10 days	

# Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
LPLFX	MP042	Specimen Type	Plain Text	Yes



MAYO CLINIC
LABORATORIES
Lymphoplasmacytic Lymphoma/Waldenstrom Macroglobulinemia, MYD88 L265P with Reflex to CXCR4, Varies

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
MP042	Specimen Type	Alphanumeric		31208-2
601511	LPLFX Reflex Result	Alphanumeric		82140-5
601510	Final Diagnosis	Alphanumeric		50398-7

LOINC and CPT codes are provided by the performing laboratory.

Sun	nlem	ental	Ren	ort.
Jub	DICIII	Ciitai	IVEN	יטועי

No

## **CPT Code Information:**

81305

#### **Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
CXCFX	CXCR4, Gene Mutation, Reflex			No	Yes (order CXLPL), Bill Only

# Reference Values:

MYD88 L265P: Mutation present or absent based on expected variant polymerase chain reaction product size for the MYD88 gene (NCBI accession NM\_002468.4).

CXCR4: Mutation present or absent in the test region c. 898-1059 (amino acids 300-353) of the CXCR4 gene (NCBI NM\_003467.2, GRCh37).