

Reporting Title: Tubular Phosp Reabsorption, Random**Performing Location:** Rochester**Specimen Requirements:**

Both serum and urine are required.

Specimen Type: Serum

Patient Preparation: Fasting

Collection Container/Tube:

Preferred: Serum gel

Acceptable: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL

Collection Instructions:

1. Centrifuge and aliquot serum into a plastic vial.
2. Label specimen as serum.

Specimen Type: Urine

Container/Tube: Plastic, 5-mL tube

Specimen Volume: 4 mL

Collection Instructions:

1. Collect a random urine specimen.
2. No preservative.
3. Label specimen as urine.

Specimen Minimum Volume:

Urine: 1 mL; Serum: See Specimen Required

Forms:

If not ordering electronically, complete, print, and send a Renal Diagnostics Test Request (T830) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	7 days	
	Refrigerated	7 days	
Urine	Refrigerated (preferred)	30 days	
	Frozen	14 days	
	Ambient	7 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
TRA	TRP	Numeric	%	In Process
GFRR	Random TmP/GFR	Numeric	mg/dL	103542-7
CRETR	Creatinine, Random, U	Numeric	mg/dL	2161-8
PHOS	Phosphorus (Inorganic), S Also used by tests: PHOS	Numeric	mg/dL	2777-1
ACREA	Creatinine, S Also used by tests: ACREA	Numeric	mg/dL	2160-0

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Components:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
RTRP	Tubular Phosp Reabsorption, Random			Yes	No
CRETR	Creatinine, Random, U			Yes	No
PHOS	Phosphorus (Inorganic), S			Yes	Yes
ACREA	Creatinine, S			Yes	Yes (Order CRTS1)

CPT Code Information:

82565
82570
84100
84105



Reference Values:

TUBULAR REABSORPTION OF PHOSPHORUS

>80%

(Although, tubular reabsorption of phosphorus levels must be interpreted in light of the prevailing plasma phosphorus and glomerular filtration rate.)

TUBULAR MAXIMUM PHOSPHORUS REABSORPTION/GLOMERULAR FILTRATION RATE (TmP/GFR)

2.6-4.4 mg/dL (0.80-1.35 mmol/L)

PHOSPHORUS (INORGANIC)

Males

1-4 years: 4.3-5.4 mg/dL

5-13 years: 3.7-5.4 mg/dL

14-15 years: 3.5-5.3 mg/dL

16-17 years: 3.1-4.7 mg/dL

> or =18 years: 2.5-4.5 mg/dL

Reference values have not been established for patients that are <12 months of age.

Females

1-7 years: 4.3-5.4 mg/dL

8-13 years: 4.0-5.2 mg/dL

14-15 years: 3.5-4.9 mg/dL

16-17 years: 3.1-4.7 mg/dL

> or =18 years: 2.5-4.5 mg/dL

Reference values have not been established for patients that are <12 months of age.

PHOSPHORUS, Random Urine

No established reference values

Random urine phosphorus may be interpreted in conjunction with serum phosphorus, using both values to calculate fractional excretion of phosphorus.

The calculation for fractional excretion (FE) of phosphorus (P) is

$$FE(P) = \frac{[P(\text{urine}) \times \text{Creat}(\text{serum})]}{[P(\text{serum}) \times \text{Creat}(\text{urine})]} \times 100$$

CREATININE Serum

Males(1)

0-11 months: 0.17-0.42 mg/dL

1-5 years: 0.19-0.49 mg/dL

6-10 years: 0.26-0.61 mg/dL

11-14 years: 0.35-0.86 mg/dL

> or =15 years: 0.74-1.35 mg/dL

Females(1)

0-11 months: 0.17-0.42 mg/dL

1-5 years: 0.19-0.49 mg/dL

6-10 years: 0.26-0.61 mg/dL

11-15 years: 0.35-0.86 mg/dL

> or =16 years: 0.59-1.04 mg/dL



CREATININE, Random Urine

16-326 mg/dL

Reference values have not been established for patients who are less than 18 years of age.