

# **Reporting Title:** Potassium, 24 HR, U **Performing Location:** Rochester

### Necessary Information:

24-Hour volume (in milliliters) is required.

### Specimen Requirements:

Supplies: Sarstedt 5 mL Aliquot Tube (T914)
Collection Container/Tube: 24-Hour graduated urine container with no metal cap or glued insert
Submission Container/Tube: Plastic, 5 mL tube or a clean, plastic aliquot container with no metal cap or glued insert
Specimen Volume: 5 mL
Collection Instructions:

Collect urine for 24 hours.
Refrigerate specimen within 4 hours of completion of 24-Hour collection.

Additional Information:
See Urine Preservatives-Collection and Transportation for 24-Hour Urine Specimens for multiple collections.

### **Specimen Minimum Volume:**

1 mL

**Toluene No** 

# **Urine Preservative Collection Options:**

Note: The addition of preservative or application of temperature controls must occur within 4 hours of completion of the collection. Ambient OK Refrigerate Preferred Frozen OK 50% Acetic Acid OK Boric Acid OK Diazolidinyl Urea OK 6M Hydrochloric Acid OK 6M Nitric Acid No Sodium Carbonate OK Thymol OK

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	30 days	

# LABORATORIES

Potassium, 24 Hour, Urine

Ambient 7 days			
Ambient 7 days			
	Ambient	7 days	
		, ,	

# Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
KUR	TM28	Collection Duration	Plain Text	Yes
KUR	VL26	Urine Volume	Plain Text	Yes

# Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
KURU	Potassium, 24 Hr, U	Numeric	mmol/24 h	2829-0
TM28	Collection Duration	Numeric	h	13362-9
VL26	Urine Volume	Numeric	mL	3167-4

LOINC and CPT codes are provided by the performing laboratory.

## Supplemental Report:

No

# **CPT Code Information:**

84133

# **Reference Values:**

> or =18 years: 16-105 mmol/24 hours

Reference values have not been established for patients who are younger than 18 years of age.