MAYO CLINIC
 Test Definition: PINZ

 LABORATORIES
 PTEN Hamartoma Tumor Syndrome, PTEN Full Gene Analysis, Varies

Reporting Title: PTEN Full Gene Analysis **Performing Location:** Rochester

Ordering Guidance:

For a comprehensive hereditary cancer panel that includes the PTEN gene, consider ordering 1 of the following tests: -BRGYP / Hereditary Breast/Gynecologic Cancer Panel, Varies -ENDCP / Hereditary Endocrine Cancer Panel, Varies -HPGLP / Hereditary Paraganglioma/Pheochromocytoma Panel, Varies -RENCP / Hereditary Renal Cancer Panel, Varies

Testing for the PTEN gene as part of a customized panel is available. For more information see CGPH / Custom Gene Panel, Hereditary, Next-Generation Sequencing, Varies.

Targeted testing for familial variants (also called site-specific or known mutations testing) is available for the genes on this panel. See FMTT / Familial Variant, Targeted Testing, Varies. To obtain more information about this testing option, call 800-533-1710.

Shipping Instructions:

Specimen preferred to arrive within 96 hours of collection.

Specimen Requirements:

Patient Preparation: A previous bone marrow transplant from an allogenic donor will interfere with testing. Call 800-533-1710 for instructions for testing patients who have received a bone marrow transplant.
Specimen Type: Whole blood
Container/Tube:
Preferred: Lavender top (EDTA) or yellow top (ACD)
Acceptable: Any anticoagulant
Specimen Volume: 3 mL
Collection Instructions:
1. Invert several times to mix blood.
2. Send whole blood specimen in original tube. Do not aliquot.

Specimen Stability Information: Ambient (preferred) 4 days/Refrigerated

Specimen Minimum Volume:

1 mL

Forms:

 New York Clients-Informed consent is required. Document on the request form or electronic order that a copy is on file. The following documents are available:
 Informed Consent for Genetic Testing (T576)
 Informed Consent for Genetic Testing-Spanish (T826)
 Molecular Genetics: Inherited Cancer Syndromes Patient Information (T519)
 MAYO CLINIC
 Test Definition: PTNZ

 LABORATORIES
 PTEN Hamartoma Tumor Syndrome, PTEN Full Gene Analysis,

3. If not ordering electronically, complete, print, and send a Oncology Test Request (T729) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
614815	Test Description	Alphanumeric		62364-5
614816	Specimen	Alphanumeric		31208-2
614817	Source	Alphanumeric		31208-2
614818	Result Summary	Alphanumeric		50397-9
614819	Result	Alphanumeric		82939-0
614820	Interpretation	Alphanumeric		69047-9
614821	Resources	Alphanumeric		99622-3
614822	Additional Information	Alphanumeric		48767-8
614823	Method	Alphanumeric		85069-3
614824	Genes Analyzed	Alphanumeric		48018-6
614825	Disclaimer	Alphanumeric		62364-5
614826	Released By	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

Supplemental

CPT Code Information:

81321



Reference Values:

An interpretive report will be provided.