

Zinc/Creatinine Ratio, Random, Urine

Reporting Title: Zinc/Creat Ratio, Random, U **Performing Location:** Rochester

Ordering Guidance:

The use of blood testing is recommended for children younger than 10 years of age. Order ZN_S / Zinc, Serum.

Specimen Requirements:

Patient Preparation: High concentrations of barium are known to interfere with most metal tests. If barium-containing contrast media has been administered, the specimen should not be collected for at least 96 hours. Supplies: Urine Tubes, 10 mL (T068)

Collection Container/Tube: Clean, plastic urine collection container with no metal cap or glued insert Submission Container/Tube: Plastic urine tube or clean, plastic aliquot container with no metal cap or glued insert Specimen Volume: 3 mL

Collection Instructions:

1. Collect a random urine specimen.

2. See Metals Analysis Specimen Collection and Transport for complete instructions.

Specimen Minimum Volume:

2 mL

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	28 days	
	Frozen	28 days	
	Ambient	14 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
615260	Zinc/Creat Ratio, U	Numeric	mcg/g Cr	13473-4
CRETR	Creatinine, Random, U	Numeric	mg/dL	2161-8

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No



Zinc/Creatinine Ratio, Random, Urine

Components:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
ZNCU	Zinc/Creat Ratio, U			Yes	No
CRETR	Creatinine, Random, U			Yes	No

CPT Code Information:

84630 82570

Reference Values:

ZINC/CREATININE: 0-17 years: Not established > or =18 years: 89-910 mcg/g creatinine

CREATITINE:

> or =18 years old: 16-326 mg/dL
Reference values have not been established for patients who are younger than 18 years of age.