

Reporting Title: HIV-1/HIV-2 RNA Detect Prenatal, P**Performing Location:** Rochester**Ordering Guidance:**

This test is indicated for qualitative detection and differentiation of HIV-1 and HIV-2 infection in pregnant individuals with indeterminate or inconclusive HIV serologic test results or who are suspected to be in the acute/early HIV infection stage (ie, "window period"). If serologic testing has not yet been performed on the patient, order either HIVSP / HIV Antigen and Antibody Prenatal Routine Screen, Plasma or HVDSP / HIV-1 and HIV-2 Antibody Confirmation and Differentiation Prenatal, Plasma, where this test will automatically be added on depending on the results of the serologic test.

Shipping Instructions:

1. Ship specimen frozen on dry ice.
2. If shipment will be delayed for more than 24 hours, freeze plasma specimen at -20 to -80 degrees C until shipment on dry ice.

Specimen Requirements:

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube: Lavender top (EDTA)

Submission Container/Tube: Plastic vial

Specimen Volume: 1.5 mL

Collection Instructions:

1. Centrifuge blood collection tube and aliquot plasma into plastic vial per collection tube manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes).
2. Freeze aliquoted plasma for shipment.

Specimen Minimum Volume:

0.8 mL

Forms:

If not ordering electronically, complete, print, and send a Microbiology Test Request (T244) with the specimen.

Specimen Type	Temperature	Time	Special Container
Plasma EDTA	Frozen (preferred)	40 days	
	Refrigerated	6 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
616344	HIV-1 RNA	Alphanumeric		25835-0
616345	HIV-2 RNA	Alphanumeric		69353-1

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

87535

87538

Reference Values:

Undetected