MAYO CLINIC Test Definition: CTNBT LABORATORIES B-Catenin Mutation Analysis, Next-Generation Sequencing, Tumor

Reporting Title: B-Catenin Mutations Analysis, Tumor **Performing Location:** Rochester

Ordering Guidance:

Multiple oncology (cancer) gene panels are available. For more information see Hematology, Oncology, and Hereditary Test Selection Guide.

Necessary Information:

A pathology report (final or preliminary), at minimum containing the following information, must accompany specimen for testing to be performed:

- 1. Patient name
- 2. Block number-must be on all blocks, slides, and paperwork (can be handwritten on the paperwork)
- 3. Tissue collection date
- 4. Source of the tissue

Specimen Requirements:

This assay requires at least 20% tumor nuclei.

-Preferred amount of tumor area with sufficient percent tumor nuclei: tissue 216 mm(2)

-Minimum amount of tumor area: tissue 36 mm(2)

-These amounts are cumulative over up to 10 unstained slides and must have adequate percent tumor nuclei.

-Tissue fixation: 10% neutral buffered formalin, not decalcified

-For specimen preparation guidance, see Tissue Requirement for Solid Tumor Next-Generation Sequencing. In this document, the sizes are given as 4 mm x 4 mm x 10 slides as preferred: approximate/equivalent to 144 mm(2) and the minimum as 3 mm x 1 mm x 10 slides: approximate/equivalent to 36 mm(2)

Preferred:

Specimen Type: Tissue block

Collection Instructions: Submit a formalin-fixed, paraffin-embedded tissue block with acceptable amount of tumor tissue.

Acceptable:

Specimen Type: Tissue slides

Slides: 1 Stained and 10 unstained

Collection Instructions: Submit 1 slide stained with hematoxylin and eosin and 10 unstained, nonbaked slides with 5-micron thick sections of the tumor tissue.

Note: The total amount of required tumor nuclei can be obtained by scraping up to 10 slides from the same block. Additional Information: Unused unstained slides will not be returned.

Specimen Type: Cytology slides (direct smears or ThinPrep)

Slides: 1 to 3 Slides

Collection Instructions: Submit 1 to 3 slides stained and coverslipped with a preferred total of 5000 nucleated cells, or a minimum of at least 3000 nucleated cells.

Note: Glass coverslips are preferred; plastic coverslips are acceptable but will result in longer turnaround times. Additional Information: Cytology slides will not be returned.

Test Definition: CTNBT LABORATORIES B-Catenin Mutation Analysis, Next-Generation Sequencing, Tumor

Specimen Minimum Volume:

See Specimen Required

Forms:

If not ordering electronically, complete, print, and send an Oncology Test Request (T729) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)		
	Refrigerated		

Result Codes:

Result ID	Reporting Name	Type Unit		LOINC®
617921	Result	Alphanumeric		82939-0
617922	Interpretation	Alphanumeric		69047-9
617923	Additional Information	Alphanumeric		48767-8
617924	Specimen	Alphanumeric		31208-2
617925	Tissue ID	Alphanumeric		80398-1
617926	Method	Alphanumeric		85069-3
617927	Disclaimer	Alphanumeric		62364-5
617928	Released By	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

88381-Microdissection, manual 81479

Test Definition: CTNBT LABORATORIES B-Catenin Mutation Analysis, Next-Generation Sequencing, Tumor

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
SLIRV	Slide Review in MG			Yes	No

Reference Values:

An interpretive report will be provided.