

Monoclonal Protein Quantitation, 24 Hour, Urine

Reporting Title: M-protein Quantitation, 24 HR, U

Performing Location: Rochester

Shipping Instructions:

Refrigerate specimen during collection and send refrigerated.

Necessary Information:

24-Hour volume (in milliliters) is required.

Specimen Requirements:

Supplies: Urine Container, 60 mL (T313)

Submission Container/Tube: Plastic, 60-mL urine bottle

Specimen Volume: 50 mL Collection Instructions: 1. Collect urine for 24 hours.

1. Collect unine for 24 hours.

2. Aliquot between 30 mL and 50 mL urine into plastic, 60-mL urine bottle.

Additional Information: See Urine Preservatives-Collection and Transportation for 24-Hour Urine Specimens for multiple collections.

Specimen Minimum Volume:

30 mL

Forms:

If not ordering electronically, complete, print, and send a Renal Diagnostics Test Request (T830) with the specimen.

Urine Preservative Collection Options:

Note: The addition of preservative or application of temperature controls must occur within 4 hours of completion of the collection.

Ambient OK

Refrigerate Preferred

Frozen OK

50% Acetic Acid No

Boric Acid No

Diazolidinyl Urea OK

6M Hydrochloric Acid No

6M Nitric Acid No

Sodium Carbonate No

Thymol OK

Toluene No

Monoclonal Protein Quantitation, 24 Hour, Urine

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------------|----------|-------------------|
| Urine | Refrigerated (preferred) | 14 days | |
| | Frozen | 5 days | |
| | Ambient | 24 hours | |

Ask at Order Entry (AOE) Questions:

| Test ID | Question ID | Description | Туре | Reportable |
|---------|-------------|---------------------|------------|------------|
| PTU3 | TM24 | Collection Duration | Plain Text | Yes |
| PTU3 | VL92 | Urine Volume | Plain Text | Yes |

Result Codes:

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|--------------------------|--------------|---------|---------|
| TP6 | Total Protein, 24 HR, U | Numeric | mg/24 h | 2889-4 |
| | Also used by tests: PTU3 | | | |
| TM24 | Collection Duration | Numeric | h | 13362-9 |
| | Also used by tests: PTU3 | | | |
| VL92 | Urine Volume | Numeric | mL | 3167-4 |
| | Also used by tests: PTU3 | | | |
| 607970 | Albumin | Alphanumeric | mg/24 h | 6941-9 |
| 607971 | Alpha-1 globulin | Alphanumeric | mg/24 h | 6794-2 |
| 607972 | Alpha-2 globulin | Alphanumeric | mg/24 h | 6795-9 |
| 607973 | Beta globulin | Alphanumeric | mg/24 h | 94714-3 |
| 607974 | Gamma globulin | Alphanumeric | mg/24 h | 94715-0 |
| 2833 | A/G Ratio | Alphanumeric | | 44294-7 |
| 21446 | M spike | Numeric | mg/24 h | 42482-0 |
| 22307 | M spike | Numeric | mg/24 h | 42482-0 |
| 21447 | Impression | Alphanumeric | | 32210-7 |

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Monoclonal Protein Quantitation, 24 Hour, Urine

Components:

| Test ID | Reporting Name | CPT Units | CPT Code | Always Performed | Orderable Separately |
|---------|-----------------------------------|-----------|----------|---------------------|-------------------------|
| PTU3 | Protein, Total, 24 HR, U | | | Yes | Yes (Order PTU) |
| PEU | Protein Electrophoresis, 24 Hr, U | | | Yes | No |

CPT Code Information:

84156 84166 0077U (if appropriate)

Reflex Tests:

| Test ID | Reporting Name | CPT Units | CPT Code | Always Performed | Orderable Separately |
|---------|------------------------------|-----------|----------|---------------------|-------------------------|
| MPTU | M-protein Mass-Fix, 24 HR, U | | | No | No |

Result Codes for Reflex Tests:

| Test ID | Result ID | Reporting Name | Туре | Unit | LOINC® |
|---------|-----------|-------------------------------------|--------------|------|--------------------|
| MPTU | 617117 | M-protein Isotype MS, 24 HR, U | Alphanumeric | | In Process |
| MPTU | 617118 | Flag M-protein Isotype MS, 24 HR, U | Alphanumeric | | No LOINC Needed |

Reference Values:

PROTEIN, TOTAL <229 mg/24 hours

Reference values have not been established for patients who are younger than 18 years of age.

ELECTROPHORESIS, PROTEIN

The following fractions, if present, will be reported as mg/24 hours:

Albumin

Alpha-1-globulin

Alpha-2-globulin

Beta-globulin



Monoclonal Protein Quantitation, 24 Hour, Urine

Gamma-globulin

MASS-FIX M-PROTEIN ISOTYPE M-protein Isotype MS: No monoclonal protein detected

Flag M-protein Isotype MS: Negative