

Monoclonal Protein Screen, 24 Hour, Urine

# **Reporting Title:** M-protein Screen, 24 HR, U **Performing Location:** Rochester

## **Shipping Instructions:**

Refrigerate specimen during collection and send refrigerated.

#### Necessary Information:

24-Hour volume (in milliliters) is required.

#### **Specimen Requirements:**

Supplies: Urine Container, 60 mL (T313)
Submission Container/Tube: Plastic, 60-mL urine bottle
Specimen Volume: 50 mL
Collection Instructions:
1. Collect urine for 24 hours.
2. Aliquot between 30 mL and 50 mL of urine into a plastic, 60-mL urine bottle.

Additional Information: See Urine Preservatives-Collection and Transportation for 24-Hour Urine Specimens for multiple collections.

#### **Specimen Minimum Volume:**

30 mL

#### Forms:

If not ordering electronically, complete, print, and send a Renal Diagnostics Test Request (T830) with the specimen.

## **Urine Preservative Collection Options:**

Note: The addition of preservative or application of temperature controls must occur within 4 hours of completion of the collection. Ambient OK Refrigerate Preferred Frozen OK 50% Acetic Acid No Boric Acid No Diazolidinyl Urea OK 6M Hydrochloric Acid No 6M Nitric Acid No Sodium Carbonate No Thymol OK Toluene No

# LABORATORIES

# **Test Definition: SMPU**

Monoclonal Protein Screen, 24 Hour, Urine

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	5 days	
	Ambient	24 hours	

# Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
PTU3	TM24	Collection Duration	Plain Text	Yes
PTU3	VL92	Urine Volume	Plain Text	Yes

# Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
617117	M-protein Isotype MS, 24 HR, U	Alphanumeric		In Process
617118	Flag M-protein Isotype MS, 24 HR, U	Alphanumeric		No LOINC Needed
TP6	Total Protein, 24 HR, U	Numeric	mg/24 h	2889-4
	Also used by tests: PTU3			
TM24	Collection Duration	Numeric	h	13362-9
	Also used by tests: PTU3			
VL92	Urine Volume	Numeric	mL	3167-4
	Also used by tests: PTU3			

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No



Monoclonal Protein Screen, 24 Hour, Urine

# **Components:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
MPTU	M-protein Mass-Fix, 24 HR, U			Yes	No
PTU3	Protein, Total, 24 HR, U			Yes	Yes

# **CPT Code Information:**

84156 84166 Electrophoresis, protein (if appropriate) 0077U

# **Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
PEU	Protein Electrophoresis, 24 Hr, U			No	No

# **Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Туре	Unit	LOINC®
PEU	607970	Albumin	Alphanumeric	mg/24 h	6941-9
PEU	607971	Alpha-1 globulin	Alphanumeric	mg/24 h	6794-2
PEU	607972	Alpha-2 globulin	Alphanumeric	mg/24 h	6795-9
PEU	607973	Beta globulin	Alphanumeric	mg/24 h	94714-3
PEU	607974	Gamma globulin	Alphanumeric	mg/24 h	94715-0
PEU	2833	A/G Ratio	Alphanumeric		44294-7
PEU	21446	M spike	Numeric	mg/24 h	42482-0
PEU	22307	M spike	Numeric	mg/24 h	42482-0
PEU	21447	Impression	Alphanumeric		32210-7



# **Test Definition: SMPU**

Monoclonal Protein Screen, 24 Hour, Urine

## **Reference Values:**

PROTEIN, TOTAL <229 mg/24 hours

Reference values have not been established for patients younger than 18 years of age.

ELECTROPHORESIS, PROTEIN The following fractions, if present, will be reported as mg/24 hours: Albumin Alpha-1-globulin Alpha-2-globulin Beta-globulin Gamma-globulin

MASS-FIX M-PROTEIN ISOTYPE M-protein Isotype MS: No monoclonal protein detected

Flag M-protein Isotype MS: Negative