

**Reporting Title:** Oxycodone w/metabolite Conf, CoC, U  
**Performing Location:** Rochester

**Specimen Requirements:**

Supplies: Chain-of-Custody Kit (T282)

Container/Tube: Chain-of-custody kit containing the specimen containers, seals, and documentation required.

Specimen Volume: 5 mL

Collection Instructions: Collect urine specimen in the container provided, seal, and submit with the associated documentation to satisfy the legal requirements for chain-of-custody testing.

**Specimen Minimum Volume:**

2 mL

**Forms:**

1. Chain of Custody Request is included in the Chain-of-Custody Kit (T282).
2. If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	28 days	
	Frozen	28 days	
	Ambient	14 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
36300	Oxycodone Immunoassay Screen	Alphanumeric	ng/mL	19642-8
61728	Oxycodone-by LC-MS/MS	Alphanumeric	ng/mL	16249-5
42012	Noroxycodone-by LC-MS/MS	Alphanumeric	ng/mL	61425-5
35964	Oxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	17395-5
42013	Noroxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	90894-7
36023	Oxycodone Interpretation	Alphanumeric		69050-3
36028	Chain of Custody	Alphanumeric		77202-0

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

80365  
G0480 (if appropriate)

**Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
COCH	Chain of Custody Processing			Yes	No
ADLTX	Adulterants Survey, CoC, U			Yes	Yes

**Result Codes for Reflex Tests:**

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
ADLTX	36121	Creatinine, U	Numeric	mg/dL	2161-8
ADLTX	36122	Specific Gravity	Numeric		5810-7
ADLTX	36123	pH	Numeric		2756-5
ADLTX	36124	Oxidants	Alphanumeric		58714-7
ADLTX	36125	Nitrites	Alphanumeric		32710-6
ADLTX	36126	Comment	Alphanumeric		48767-8

**Reference Values:**

Negative

Cutoff concentrations:

Oxycodone Immunoassay screen: 100 ng/mL

By liquid chromatography-tandem mass spectroscopy:

Oxycodone: 25 ng/mL

Noroxycodone: 25 ng/mL

Oxymorphone: 25 ng/mL

Noroxymorphone: 25 ng/mL