

Powassan Virus, RNA, Molecular Detection, PCR, Serum

Reporting Title: Powassan Virus, PCR, Serum **Performing Location:** Rochester

Ordering Guidance:

The Powassan virus polymerase chain reaction result may remain positive for a longer time in urine than in blood, serum, and cerebrospinal fluid (7 days or more).

Additional Testing Requirements:

This test is designed to be used during the first 7 days after symptom onset and should be ordered in conjunction with POWV / Powassan Virus, IgM, Enzyme-Linked Immunosorbent Assay, Serum.

Specimen Requirements:

Collection Container/Tube: Preferred: Serum gel Acceptable: Red top Submission Container: Plastic vial Specimen Volume: 0.5 mL Collection Instructions: 1. Centrifuge and aliquot serum into a plastic vial within 2 hours of collection. 2. Label specimen as serum

Specimen Minimum Volume:

0.3 mL

Forms:

If not ordering electronically, complete, print, and send a Microbiology Test Request (T732) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	7 days	



Powassan Virus, RNA, Molecular Detection, PCR, Serum

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
617490	Lineage 1	Alphanumeric		55144-0
618146	Lineage 2 (Deer Tick Virus)	Alphanumeric		55144-0
618150	Powassan Virus PCR Comment	Alphanumeric		77202-0

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

87798

Reference Values:

Negative