

**Reporting Title:** Orthostatic Proteinuria, Random, U**Performing Location:** Rochester**Specimen Requirements:**

Patient Preparation: Specimens should be collected before fluorescein is given or not collected until at least 24 hour later.

Container/Tube: Plastic, 5-mL tube

Specimen Volume: 4 mL

Collection Instructions:

1. This is a nighttime (supine) collection.
2. Void immediately prior to lying down to sleep.
3. Collect an 8-hour (nighttime) urine specimen.
4. Upon awaking in the morning, the patient should immediately void into the specimen cup.
4. No preservative.
5. Invert well before taking 4-mL (maximum) aliquot.

**Specimen Minimum Volume:**

1 mL

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	30 days	
	Ambient	24 hours	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
PRCON	Protein, Total, Random, U	Numeric	mg/dL	2888-6
RATO2	Protein/Creatinine Ratio	Numeric	mg/mg	2890-2
CRETR	Creatinine, Random, U	Numeric	mg/dL	2161-8

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

---

**Components:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
PRCON	Protein, Total, Random, U			Yes	No
RATO2	Protein/Creatinine Ratio			Yes	No
CRETR	Creatinine, Random, U			Yes	No

**CPT Code Information:**

CRETR- 82570  
PRCON- 84156

**Reference Values:**

PROTEIN/CREATININE RATIO:

> or =18 years: <0.18 mg/mg creatinine

CREATININE:> or =18 years = 16-326 mg/dL

Reference values have not been established for patients younger than 18 years of age.