

## **Test Definition: VZIKS**

Zika Virus, PCR, Molecular Detection, Serum

Reporting Title: Zika Virus PCR, Serum

Performing Location: Rochester

#### Ordering Guidance:

Due to similar clinical presentations, testing for RNA or IgM-class antibodies to dengue and chikungunya viruses, concurrently with Zika virus testing, is strongly recommended.

## **Specimen Requirements:**

Collection Container/Tube: Preferred: Serum gel Acceptable: Red top

Submission Container: Sterile container

Specimen Volume: 0.5 mL Collection Instructions:

1. Collect whole blood in a serum gel tube.

2. Centrifuge and aliquot the serum into a sterile container within 6 hours of collection.

3. Label specimen as serum.

#### **Specimen Minimum Volume:**

0.25 mL

#### Forms:

If not ordering electronically, complete, print, and send a Microbiology Test Request (T244) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	7 days	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
619456	Zika Serum PCR Result	Alphanumeric		85622-9
619457	Zika Serum PCR Interpretation	Alphanumeric		69048-7

LOINC and CPT codes are provided by the performing laboratory.

## Supplemental Report:



# **Test Definition: VZIKS**

Zika Virus, PCR, Molecular Detection, Serum

N I	_
IN	U

## **CPT Code Information:**

87662

## **Reference Values:**

Negative

Reference values apply to all ages.