

Reporting Title: HPV with Genotyping, PCR, ThinPrep**Performing Location:** Rochester**Ordering Guidance:**

For vaginal specimens, order VHPV / Human Papillomavirus (HPV) Vaginal Detection with Genotyping for High-Risk Types by PCR.

Necessary Information:

Specimen source is required, collection date, and patient identifiers are required.

Specimen Requirements:

Specimen Type: Cervical (endocervical or ectocervical)

Container/Tube: ThinPrep/PreservCyt solution vial

Specimen Volume: 3 mL of solution in ThinPrep/PreservCyt vial

Collection Instructions:

1. Bag ThinPrep specimens individually as they have a tendency to leak during transport.
2. Place labels on the vial and on the bag.

Specimen Minimum Volume:

1 mL

Forms:

If not ordering electronically, complete, print, and send Microbiology Test Request (T244) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)	42 days	
	Refrigerated	42 days	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Reportable
---------	-------------	-------------	------	------------

Test ID	Question ID	Description	Type	Reportable
HPV	SS017	Specimen Source: <ul style="list-style-type: none">• Cervix• Cervix Biopsy• Endocervix• Cervix / Endocervix	Answer List	Yes

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
SS017	Specimen Source	Alphanumeric		31208-2
35924	HPV High Risk type 16, PCR	Alphanumeric		77399-4
35925	HPV High Risk type 18, PCR	Alphanumeric		77400-0
35926	HPV other High Risk types, PCR	Alphanumeric		71431-1

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

87624

G0476 (if appropriate)

Reference Values:

Negative for human papillomavirus (HPV) genotypes 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68