Test Definition: AMPMX

Amphetamine-Type Stimulants Confirmation, Chain of Custody, Meconium

Reporting Title: Amphetamines Confirmation, CoC, M

Performing Location: Rochester

Specimen Requirements:

Supplies: Chain of Custody Meconium Kit (T653) includes the specimen containers, seals, and documentation required.

Specimen Volume: 1 g (approximately 1 teaspoon)

Collection Instructions: Collect entire random meconium specimen.

Additional Information:

1. Specimen that arrives with a broken seal does not meet the chain of custody requirements.

2. The laboratory recommends sending chain-of-custody specimens by overnight shipment.

Specimen Minimum Volume:

0.3 g (approximately 1/4 teaspoon)

Forms:

- 1. Chain of Custody Request is included in the Chain-of-Custody Meconium Kit (T653).
- 2. If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Meconium	Frozen (preferred)	28 days	
	Ambient	28 days	
	Refrigerated	28 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
36136	Amphetamine	Alphanumeric	ng/g	26959-7
36137	Methamphetamine	Alphanumeric	ng/g	69022-2
36138	3,4-methylenedioxyamphetamine	Alphanumeric	ng/g	69023-0
36139	3,4-methylenedioxyethylamphetamine	Alphanumeric	ng/g	69024-8
36140	3,4-methylenedioxymethamphetamine	Alphanumeric	ng/g	69025-5
36141	Interpretation	Alphanumeric		69050-3
36142	Chain of Custody	Alphanumeric		77202-0

LOINC and CPT codes are provided by the performing laboratory.



Test Definition: AMPMX

Amphetamine-Type Stimulants Confirmation, Chain of Custody, Meconium

Supplemental Report:

No

CPT Code Information:

80324 80359 G0480 (if appropriate)

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
СОСН	Chain of Custody Processing			Yes	No

Reference Values:

Negative

Positives are reported with a quantitative liquid chromatography tandem mass spectrometry (LC-MS/MS) result. Cutoff concentrations for LC-MS/MS testing:

Amphetamine: 20 ng/g Methamphetamine: 20 ng/g

3,4-Methylenedioxyamphetamine: 20 ng/g 3,4-Methylenedioxyethylamphetamine: 20 ng/g 3,4-Methylenedioxymethamphetamine: 20 ng/g