

## **Test Definition: NGSHM**

MayoComplete Myeloid Neoplasms, Comprehensive OncoHeme Next-Generation Sequencing, Varies

Reporting Title: Myeloid Neoplasms, NGS, V

Performing Location: Rochester

## **Shipping Instructions:**

Peripheral blood and bone marrow specimens must arrive within 14 days of collection.

## **Necessary Information:**

The following information is required:

- 1. Clinical diagnosis
- Pertinent clinical history, including disease phase (diagnostic, remission, relapse/refractory) and therapy status (especially if patient has received a hematopoietic stem cell transplant).
- 3. Clinical or morphologic suspicion
- 4. Date of collection
- 5. Specimen source

## **Specimen Requirements:**

Submit only 1 of the following specimens:

Preferred Specimen Type: Bone marrow aspirate

Container/Tube:

Preferred: Lavender top (EDTA) or yellow top (ACD)

Acceptable: Green top (sodium heparin)

Specimen Volume: 2 mL Collection Instructions:

- 1. Invert several times to mix bone marrow.
- 2. Send bone marrow specimen in original tube. Do not aliquot.
- 3. Label specimen as bone marrow.

Specimen Stability: Ambient (preferred)/Refrigerate

Specimen Type: Whole blood

Container/Tube:

Preferred: Lavender top (EDTA) or yellow top (ACD)

Acceptable: Green top (sodium heparin)

Specimen Volume: 3 mL Collection Instructions:

- 1. Invert several times to mix blood.
- 2. Send whole blood specimen in original tube. Do not aliquot.
- 3. Label specimen as blood.

Specimen Stability: Ambient (preferred)/Refrigerate

Specimen Type: Extracted DNA from blood or bone marrow

Container/Tube: 1.5- to 2-mL tube with indication of volume and concentration of the DNA

Specimen Volume: Entire specimen

Collection Instructions: Label specimen as extracted DNA and source of specimen

Specimen Stability: Frozen (preferred)/Refrigerate/Ambient



# **Test Definition: NGSHM**

MayoComplete Myeloid Neoplasms, Comprehensive OncoHeme Next-Generation Sequencing, Varies

## **Specimen Minimum Volume:**

Blood, Bone marrow: 1 mL

Extracted DNA: 100 mcL at 20 ng/mcL concentration

#### Forms:

1. Hematopathology Patient Information (T676)

2. If not ordering electronically, complete, print, and send a Hematopathology/Cytogenetics Test Request (T726) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies	14 days	

## Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
NGSHM	MP024	Specimen Type	Plain Text	Yes
NGSHM	NGSD	Indication for Test	Plain Text	Yes

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
MP024	MP024 Specimen Type			31208-2
NGSD	Indication for Test	Alphanumeric		42349-1
601696	NGSHM Result	Alphanumeric		No LOINC Needed
37276	Pathogenic Mutations Detected	Alphanumeric		82939-0
37283	Interpretation	Alphanumeric		69047-9
37282	Clinical Trials	Alphanumeric		82786-5
37277	Variants of Unknown Significance	Alphanumeric		93367-1
37278	Additional Notes	Alphanumeric		48767-8
37279	Method Summary	Alphanumeric		85069-3
37420	Disclaimer	Alphanumeric		62364-5



# **Test Definition: NGSHM**

MayoComplete Myeloid Neoplasms, Comprehensive OncoHeme Next-Generation Sequencing, Varies

Result ID	Reporting Name	Туре	Unit	LOINC®
37280	OncoHeme Panel Gene list	Alphanumeric		36908-2
37287	Reviewed By:	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

<b>Supplemental</b>	Report:
No	

## **CPT Code Information:**

81450

#### **Reference Values:**

An interpretive report will be provided.