

---

**Reporting Title:** NR4A3, Rearrangement, FISH, Ts  
**Performing Location:** Rochester

**Shipping Instructions:**

Advise Express Mail or equivalent if not on courier service.

**Necessary Information:**

1. A pathology report is required in order for testing to be performed. Acceptable pathology reports include working drafts, preliminary pathology or surgical pathology reports.
2. A reason for testing must be provided. If this information is not provided, an appropriate indication for testing may be entered by Mayo Clinic Laboratories.

**Specimen Requirements:**

Submit only 1 of the following specimens:

Specimen Type: Tissue

Container/Tube: Formalin-fixed, paraffin-embedded tumor tissue block

Specimen Type: Slides

Specimen Volume: Four Consecutive, unstained, 5 micron-thick sections placed on positively charged slides and 1 hematoxylin and eosin-stained slide

**Specimen Minimum Volume:**

Two consecutive, unstained, 5 micron-thick sections placed on positively charged slides and 1 hematoxylin and eosin-stained slide.

Specimen Type	Temperature	Time	Special Container
Tissue	Ambient (preferred)		
	Refrigerated		

**Ask at Order Entry (AOE) Questions:**

Test ID	Question ID	Description	Type	Reportable
NR4A3	CG999	Reason for Referral	Plain Text	Yes

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
92340	Result Summary	Alphanumeric		50397-9
92341	Interpretation	Alphanumeric		69965-2
92342	Result	Alphanumeric		62356-1
CG999	Reason For Referral	Alphanumeric		42349-1
92343	Specimen	Alphanumeric		31208-2
92344	Source	Alphanumeric		31208-2
92345	Tissue ID	Alphanumeric		80398-1
92346	Method	Alphanumeric		85069-3
92347	Additional Information	Alphanumeric		48767-8
92348	Disclaimer	Alphanumeric		62364-5
92349	Released By	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

88271 x 2, 88291-DNA probe, each (first probe set), Interpretation and report  
 88271 x 2-NA probe, each; each additional probe set (if appropriate)  
 88271 x 1-DNA probe, each; coverage for sets containing 3 probes (if appropriate)  
 88271 x 2-DNA probe, each; coverage for sets containing 4 probes (if appropriate)  
 88271 x 3-DNA probe, each; coverage for sets containing 5 probes (if appropriate)  
 88274 w/modifier 52-Interphase in situ hybridization, &lt;25 cells, each probe set (if appropriate)  
 88274-Interphase in situ hybridization, 25 to 99 cells, each probe set (if appropriate)  
 88275-Interphase in situ hybridization, 100 to 300 cells, each probe set (if appropriate)

**Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
_PBCT	Probe, +2			No	No (Bill Only)

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
_PADD	Probe, +1			No	No (Bill Only)
_PB02	Probe, +2			No	No (Bill Only)
_PB03	Probe, +3			No	No (Bill Only)
_IL25	Interphases,			No	No (Bill Only)
_I099	Interphases, 25-99			No	No (Bill Only)
_I300	Interphases, >=100			No	No (Bill Only)

### Reference Values:

An interpretive report will be provided.