

## **Test Definition: CZPS**

Clonazepam and 7-Aminoclonazepam, Serum

Reporting Title: Clonazepam and 7-Aminoclonazepam, S

Performing Location: Rochester

## **Specimen Requirements:**

Collection Container/Tube: Red top Submission Container/Tube: Plastic vial

Specimen Volume: 1.2 mL Collection Instructions:

1. Draw blood immediately before next scheduled dose (minimum 12 hours after last dose).

2. Within 2 hours of collection, the specimen must be centrifuged and the serum aliquoted into a plastic vial.

## **Specimen Minimum Volume:**

0.6 mL

#### Forms:

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container	
Serum Red	Refrigerated (preferred)	14 days		
	Frozen	28 days		
	Ambient	72 hours		

#### Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
65044	Clonazepam	Numeric	ng/mL	3494-2
41782	7-Aminoclonazepam	Numeric	ng/mL	28059-4

LOINC and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

## **CPT Code Information:**

80346

G0480 (if appropriate)



# **Test Definition: CZPS**

Clonazepam and 7-Aminoclonazepam, Serum

## **Reference Values:**

Clonazepam

Anticonvulsant: 20-70 ng/mL Anxiolytic: 4-80 ng/mL

Some individuals may show therapeutic response outside of these ranges or may display toxicity within the therapeutic range, thus interpretation should include clinical evaluation.

Note: Therapeutic ranges are for specimens collected at trough (ie, immediately before next scheduled dose). Levels may be elevated in non-trough specimens.