

Test Definition: KKRP

Kingella kingae, Molecular Detection, PCR, Varies

Reporting Title: Kingella kingae PCR

Performing Location: Rochester

Necessary Information:

Specimen source is required.

Specimen Requirements:

The high sensitivity of amplification by polymerase chain reaction requires the specimen to be processed in an environment in which contamination of the specimen by Kingella kingae DNA is unlikely.

Submit only 1 of the following specimens:

Specimen Type: Synovial fluid Preferred: Lavender top (EDTA)

Acceptable: Pink top (EDTA), royal blue top (EDTA), sterile vial containing EDTA-derived aliquot, red clot tube (no

anticoagulant), or sterile container

Specimen Volume: 0.5 mL

Collection Instructions: Send specimen in original tube (preferred).

Specimen Stability Information: Refrigerated (preferred) <7 days /Frozen <7 days

Specimen Type: Fresh tissue or biopsy

Sources: Bone, joint, synovium, heart valve, aorta, or endocardium

Container/Tube: Sterile container

Specimen Volume: Entire collection or 5 mm(3)- approximately the size of a pencil eraser

Collection Instructions:

- 1. Collect fresh tissue specimen.
- Submit tissue only, do not add fluid to tissue
- 3. Refrigerate or freeze specimen.

Specimen Stability Information: Refrigerated (preferred) <7 days/ Frozen <7 days

Preferred Paraffin-embedded tissue block:

Specimen Type: Formalin-fixed, paraffin-embedded tissue block (FFPE) Sources: Bone, joint, synovium, heart valve, aorta, or endocardium

Supplies: Tissue Block Container (T553)

Container/Tube: Tissue block

Collection Instructions: Submit a formalin-fixed, paraffin-embedded tissue block to be cut and returned.

Specimen Stability Information: Ambient (preferred)/Refrigerated

Acceptable Paraffin-embedded tissue block:

Specimen Type: Formalin-fixed, paraffin-embedded tissue block (FFPE) Sources: Bone, joint, synovium, heart valve, aorta, or endocardium Container/Tube: Sterile container for each individual cut section (scroll).

Collection Instructions: Perform microtomy and prepare five separate 10-micron sections. Each section (scroll) must be

placed in a separate sterile container for submission.

Specimen Stability Information: Ambient (preferred)/Refrigerated



Test Definition: KKRP

Kingella kingae, Molecular Detection, PCR, Varies

Specimen Minimum Volume:

Fluid/fresh tissue or biopsy: See Specimen Required Paraffin-embedded tissue block: Two 10-micron sections

Forms:

If not ordering electronically, complete, print, and send a Microbiology Test Request (T244) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
KKRP	KKSRC	Specimen Source	Plain Text	Yes

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
KKSRC	Specimen Source	Alphanumeric		31208-2
48324	Kingella kingae PCR	Alphanumeric		65809-6

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

87798



Test Definition: KKRP

Kingella kingae, Molecular Detection, PCR, Varies

R	efe	ren	ce V	/alı	ues:

Not applicable