

# **Test Definition: HIQNP**

HIV-1 RNA Detection and Quantification, Prenatal, Plasma

Reporting Title: HIV-1 RNA Detect/Quant Prenatal, P

Performing Location: Rochester

## Ordering Guidance:

For detection and quantification of HIV-1 RNA in plasma before and during antiretroviral therapy in pregnant individuals only. For nonpregnant individuals, order HIVQN / HIV-1 RNA Detection and Quantification, Plasma.

### **Shipping Instructions:**

Ship specimen frozen on dry ice. If shipment will be delayed for more than 24 hours, freeze plasma specimen at -20 to -80 degrees C until shipment

## Specimen Requirements:

Supplies: Sarstedt Aliquot Tube, 5 mL (T914) Collection Container/Tube: Lavender top (EDTA)

Submission Container/Tube: Plastic vial

Specimen Volume: 1.5 mL Collection Instructions:

- 1. Centrifuge blood collection tube and aliquot plasma into plastic vial per collection tube manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes).
- 2. Freeze aliquoted plasma for shipment.

## **Specimen Minimum Volume:**

0.8 mL

#### Forms:

If not ordering electronically, complete, print, and send a Microbiology Test Request (T732) with the specimen.

Specimen Type	Temperature	Time	Special Container
Plasma EDTA	Frozen (preferred)	84 days	
	Refrigerated	6 days	_



# **Test Definition: HIQNP**

HIV-1 RNA Detection and Quantification, Prenatal, Plasma

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
65567	HIV-1 RNA Detect/Quant Prenatal, P	Alphanumeric	copies/mL	70241-5

LOINC and CPT codes are provided by the performing laboratory.

Supp	lemental	Report:
------	----------	---------

No

# **CPT Code Information:**

87536

# **Reference Values:**

Undetected