MAYO CLINIC
 Test Definition: FLLAG

 LABORATORIES
 Certolizumab and Anti-Certolizumab Antibody, DoseASSURE CTZ

## **Reporting Title:** Certolizumab and Anti-Certo Ab Performing Location: Esoterix Endocrinology

#### Specimen Requirements:

Specimen Type: Serum Container/Tube: Red or SST Specimen Volume: 2 mL Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 2 mL of serum frozen in a plastic vial. To avoid delays in turnaround time when requesting multiple tests, please submit separate frozen specimens for each test requested.

#### Specimen Minimum Volume:

0.60 mL (Note: This volume does not allow for repeat testing.)

#### Forms:

If not ordering electronically, complete, print, and send Gastroenterology and Hepatology Test Request (T728) with the specimen

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	14 days	
	Refrigerated	14 days	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z5637	Certolizumab	Alphanumeric		87404-0
Z5638	Anti-Certolizumab Antibody	Alphanumeric		87405-7

LOINC and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

#### **CPT Code Information:**

80299 82397

# Test Definition: FCZAC LABORATORIES Certolizumab and Anti-Certolizumab Antibody, DoseASSURE CTZ

### **Reference Values:**

Certolizumab: Quantitation Limit: <1.0 ug/mL

Results of 1 ug/mL or higher indicate detection of certolizumab

Anti-Certolizumab Antibody: Quantitation Limit: <40 ng/mL

Results of 40 ng/mL or higher indicate detection of anti-certolizumab pegol antibodies.