

Reporting Title: Fluoxetine, S**Performing Location:** Rochester**Specimen Requirements:**

Collection Container/Tube: Red top (serum gel/SST are not acceptable)

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions:

1. Draw blood immediately before the next scheduled dose (trough).
2. Centrifuge and aliquot serum into plastic vial within 2 hours of collection.

Specimen Minimum Volume:

0.5 mL

Forms:

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	28 days	
	Frozen	28 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
80228	Fluoxetine, S	Numeric	ng/mL	74982-0
251	Norfluoxetine, S	Numeric	ng/mL	3868-7
252	Fluoxetine+Norfluoxetine	Alphanumeric	ng/mL	74948-1

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80299

Reference Values:

Fluoxetine + Norfluoxetine: 120-500 ng/mL