

Barbiturates Confirmation, Random, Urine

# **Reporting Title:** Barbiturates Confirmation, U **Performing Location:** Rochester

### Ordering Guidance:

For situations where chain of custody is required, a Chain-of-Custody Kit (T282) is available. For chain-of-custody testing, order BARBX / Barbiturates Confirmation, Chain of Custody, Random, Urine.

Additional drug panels and specific requests are available. Call 800-533-1710 or 507-266-5700.

### Additional Testing Requirements:

If urine creatinine is required or adulteration of the sample is suspected, also order ADULT / Adulterants Survey, Random, Urine in addition to this test.

### **Specimen Requirements:**

Supplies: Urine Tubes, 10 mL (T068)
Container/Tube: Plastic, 10-mL urine tube
Specimen Volume: 10 mL
Collection Instructions:
1. Collect a random urine specimen.
2. No preservative.
Additional Information:
1. No specimen substitutions.
2. STAT requests are not accepted for this test.

### **Specimen Minimum Volume:**

5 mL

### Forms:

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	14 days	
	Ambient	72 hours	



# **Test Definition: BARBU**

Barbiturates Confirmation, Random, Urine

# **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
21442	Butalbital-by GC/MS	Alphanumeric	ng/mL	16237-0
20642	Amobarbital-by GC/MS	Alphanumeric	ng/mL	16239-6
20643	Pentobarbital-by GC/MS	Alphanumeric	ng/mL	16240-4
20644	Secobarbital-by GC/MS	Alphanumeric	ng/mL	16238-8
20645	Phenobarbital-by GC/MS	Alphanumeric	ng/mL	16241-2
20651	Barbiturates Interpretation	Alphanumeric		69050-3
20647	Chain of Custody	Alphanumeric		77202-0

LOINC and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

# **CPT Code Information:**

80345 G0480 (if appropriate)

## **Reference Values:**

Negative

Cutoff concentrations: BUTALBITAL BY GC-MS <100 ng/mL

AMOBARBITAL BY GC-MS <100 ng/mL

PENTOBARBITAL BY GC-MS <100 ng/mL

SECOBARBITAL BY GC-MS <100 ng/mL

PHENOBARBITAL BY GC-MS <100 ng/mL