

Proinsulin, Plasma

# **Reporting Title:** Proinsulin, P **Performing Location:** Rochester

### **Specimen Requirements:**

Patient Preparation:
1.Patient should be fasting for 8 hours.
2. Infants younger than 2 years of age should fast a maximum of 6 hours.
Supplies: Sarstedt Aliquot Tube, 5 mL (T914)
Collection Container/Tube: Ice-cooled, lavender top (EDTA)
Submission Container/Tube: Plastic vial
Specimen Volume: 0.5 mL
Collection Instructions:
1. After collection, place the whole blood on ice for at least 10 minutes, then centrifuge at refrigerated temperature.
2. Aliquot plasma into a plastic vial and send frozen.

### **Specimen Minimum Volume:**

0.25 mL

Specimen Type	Temperature	Time	Special Container	
Plasma EDTA	Frozen	30 days		

## Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
80908	Proinsulin, P	Numeric	pmol/L	27882-0

LOINC and CPT codes are provided by the performing laboratory.

## Supplemental Report:

No

# **CPT Code Information:**

84206



# **Test Definition: PINS**

Proinsulin, Plasma

#### **Reference Values:**

3.6-22 pmol/L