

Reporting Title: Digoxin, Free, S**Performing Location:** Rochester**Specimen Requirements:**

Patient Preparation: For 12 hours before specimen collection do not take multivitamins or dietary supplements containing biotin (vitamin B7), which is commonly found in hair, skin, and nail supplements and multivitamins.

Collection Container/Tube:

Preferred: Serum gel

Acceptable: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions:

1. Draw blood 6 to 8 hours after last dose of digoxin.
2. Serum gel tubes should be centrifuged within 2 hours of collection.
3. Red-top tubes should be centrifuged, and the serum aliquoted into a plastic vial within 2 hours of collection.

Specimen Minimum Volume:

0.6 mL

Forms:

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

-Cardiovascular Test Request Form (T724)

-Therapeutics Test Request (T831)

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	180 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FRDIG	Digoxin, Free, S	Numeric	ng/mL	3562-6

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80163

Reference Values:

<16 years:

Therapeutic ranges have not been established for patients who are under 16 years of age. In adults, the suggested serum free digoxin therapeutic range is 0.4-0.9 ng/mL.

Toxic concentration: > or =3.0

> or =16 years:

0.4-0.9 ng/mL

Toxic concentration: > or =3.0 ng/ mL