

Test Definition: PKU

Phenylalanine and Tyrosine, Plasma

Reporting Title: Phenylalanine and Tyrosine, P

Performing Location: Rochester

Necessary Information:

1. Patient's age is required.

2. Include family history, clinical condition (asymptomatic or acute episode), diet, and drug therapy information.

Specimen Requirements:

Patient Preparation: Fasting (4 hours or more for infants)

Collection Container/Tube:

Preferred: Green top (sodium heparin)

Acceptable: Green top (lithium heparin), Lavender top (EDTA)

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL Collection Instructions:

1. Centrifuge specimen and aliquot plasma into plastic vial.

2. Send plasma frozen.

Specimen Minimum Volume:

0.1 mL

Forms:

If not ordering electronically, complete, print, and send a Biochemical Genetics Test Request (T798) with the specimen.

Specimen Type	Temperature	Time	Special Container
Plasma	Frozen (preferred)	14 days	
	Refrigerated	14 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
8380	Phenylalanine, P	Numeric	nmol/mL	14875-9
8627	Tyrosine, P	Numeric	nmol/mL	20660-7

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:



Test Definition: PKU

Phenylalanine and Tyrosine, Plasma

No

CPT Code Information:

84030-Phenylalanine 84510-Tyrosine 82542 (if appropriate for government payers)

Reference Values:

PHENYLALANINE

Premature: 98-213 nmol/mL 0-31 days: 38-137 nmol/mL 1-24 months: 31-75 nmol/mL 2-18 years: 26-91 nmol/mL > or =19 years: 35-85 nmol/mL

Conversion Formulas:

Result in mg/dL x 60.5=result in nmol/mL Result in nmol/mL x 0.0165=result in mg/dL

TYROSINE

Premature: 147-420 nmol/mL 0-31 days: 55-147 nmol/mL 1-24 months: 22-108 nmol/mL 2-18 years: 24-115 nmol/mL > or =19 years: 34-112 nmol/mL

Conversion Formulas:

Result in mg/dL x 55.2=result in nmol/mL Result in nmol/mL x 0.0181=result in mg/dL