

Opiate Confirmation, Meconium

# **Reporting Title:** Opiate Confirmation, M **Performing Location:** Rochester

### **Ordering Guidance:**

For chain-of-custody testing, order OPTMX / Opiate Confirmation, Chain of Custody, Meconium.

### Specimen Requirements:

Supplies: Stool container, Small (Random), 4 oz (T288) Container/Tube: Stool container Specimen Volume: 1 g (approximately 1 teaspoon) Collection Instructions: Collect entire random meconium specimen.

## Specimen Minimum Volume:

0.3 g (approximately 1/4 teaspoon)

#### Forms:

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Meconium	Frozen (preferred)	28 days	
	Refrigerated	28 days	
	Ambient	14 days	

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
31848	Morphine	Alphanumeric	ng/g	69027-1
31850	Oxymorphone	Alphanumeric	ng/g	69028-9
31849	Hydromorphone	Alphanumeric	ng/g	68541-2
31847	Codeine	Alphanumeric	ng/g	68542-0
31852	Oxycodone	Alphanumeric	ng/g	68543-8
31851	Hydrocodone	Alphanumeric	ng/g	68544-6
31868	Interpretation	Alphanumeric		8215-6



**Opiate Confirmation, Meconium** 

Result ID	Reporting Name	Туре	Unit	LOINC®
31869	Chain of Custody	Alphanumeric		77202-0

LOINC and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

# **CPT Code Information:**

80361 80365 G0480 (if appropriate)

# **Reference Values:**

Negative

Positives are reported with a quantitative liquid chromatography tandem mass spectrometry (LC-MS/MS) result. Cutoff concentrations for LC-MS/MS testing: Codeine: 20 ng/g Hydrocodone: 20 ng/g Hydromorphone: 20 ng/g Morphine: 20 ng/g Oxycodone: 20 ng/g Oxymorphone: 20 ng/g