

**Test Definition: THEO** 

Theophylline, Serum

**Reporting Title:** Theophylline, S **Performing Location:** Rochester

# **Specimen Requirements:**

Container/Tube: Preferred: Serum gel Acceptable: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL Collection Instructions:

1. Serum gel tubes should be centrifuged within 2 hours of collection.

2. Red-top tubes should be centrifuged and the serum aliquoted into a plastic vial within 2 hours of collection.

## **Specimen Minimum Volume:**

0.25 mL

### Forms:

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Serum         | Refrigerated (preferred) | 14 days  |                   |
|               | Frozen                   | 28 days  |                   |
|               | Ambient                  | 72 hours |                   |

### **Result Codes:**

| Result ID | Reporting Name  | Туре    | Unit   | LOINC® |
|-----------|-----------------|---------|--------|--------|
| 8661      | Theophylline, S | Numeric | mcg/mL | 4049-3 |

LOINC and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

### **CPT Code Information:**

80198



**Test Definition: THEO** 

Theophylline, Serum

# **Reference Values:**

Therapeutic:

Bronchodilation: 8.0-20.0 mcg/mL

Neonatal apnea (< or =4 weeks old): 6.0-13.0 mcg/mL Critical value: >20.0 mcg/mL