
Reporting Title: JAK2 Exon 12 Mutation Detection, B**Performing Location:** Rochester**Ordering Guidance:**

In all cases being evaluated for JAK2 mutation status, the initial test that should be ordered is JAK2B / JAK2 V617F Mutation Detection, Blood, a sensitive assay for detection of the mutation. However, if no JAK2 V617F mutation is found, further evaluation of JAK2 may be clinically indicated.

Shipping Instructions:

Specimen must arrive within 5 days (120 hours) of collection. Collect and package specimen as close to shipping time as possible.

Necessary Information:

Date of collection is required.

Specimen Requirements:

Container/Tube:

Preferred: Lavender top (EDTA)

Acceptable: Yellow top (ACD)

Specimen Volume: 10 mL

Collection Instructions:

1. Invert several times to mix blood.
2. Send whole blood specimen in original tube. Do not aliquot.

Specimen Minimum Volume:

8 mL

Forms:

1. Hematopathology Patient Information (T676)
2. If not ordering electronically, complete, print, and send a Hematopathology/Cytogenetics Test Request (T726) with the specimen.

Specimen Type	Temperature	Time	Special Container
Whole blood	Refrigerated (preferred)	5 days	PURPLE OR PINK TOP/EDTA

	Ambient	5 days	PURPLE OR PINK TOP/EDTA
--	---------	--------	----------------------------

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
39467	JAK2 Sequencing Result	Alphanumeric		55300-8
20194	Final Diagnosis:	Alphanumeric		34574-4

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

0027U-JAK2 (Janus kinase 2) (eg, myeloproliferative disorder), exon 12 sequence and exon 13 sequence, if performed

Reference Values:

An interpretive report will be provided.