

---

**Reporting Title:** Phenosense Entry Fuzeon  
**Performing Location:** Monogram Biosciences, Inc**Specimen Requirements:**

Specimen Type: Plasma

Container/Tube: Lavender-top (EDTA) or plasma preparation tube (PPT)

Specimen Volume: 3 mL

Collection Instructions: Draw blood into two 5-mL PPT (pearl top) or EDTA (lavender top) tube. Centrifuge specimen within six hours of collection. Transfer plasma to one or more polypropylene screw-capped tube(s) and freeze. Freeze immediately. Send 3 mL plasma in a screw-cap vial frozen.

To avoid delays in turnaround time when requesting multiple tests, please submit separate frozen specimens for each test requested.

**RECOMMENDED:**

1. Patient's most recent viral load
2. Viral load collection date

**NOTE:**

1. Intended for use only for patients with viral loads greater than or equal to 500 copies/mL. For best results, viral loads should be confirmed within two weeks prior to submission for testing at Monogram.
2. Patient samples submitted <30 days apart are considered duplicate and will be canceled.

**Specimen Minimum Volume:**

1.0 mL

Specimen Type	Temperature	Time	Special Container
Plasma EDTA	Frozen	14 days	

**Ask at Order Entry (AOE) Questions:**

Test ID	Question ID	Description	Type	Reportable
FPFUZ	Z2098	Most Recent Viral Load	Plain Text	No
FPFUZ	Z2099	Viral Load Collected Date	Plain Text	No

---

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
Z1041	Phenosense Entry Fuzeon	Alphanumeric		45182-3

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**  
Referral**CPT Code Information:**

87903

**Reference Values:**

A final report will be attached in MayoAccess.