

Test Definition: MMLNS

Antimicrobial Susceptibility, Nocardia species and other Aerobic Actinomycetes, Varies

Reporting Title: Susc, Aerobic Actinomycetes

Performing Location: Rochester

Shipping Instructions:

1. For shipping information see <u>Infectious Specimen Shipping Guidelines</u>.

2. Place specimen in a large infectious container and label as an etiologic agent/infectious substance.

Necessary Information:

Specimen source is required.

Specimen Requirements:

Specimen Type: Organism in pure culture **Supplies:** Infectious Container, Large (T146)

Container/Tube:

Preferred: Middlebrook 7H10 agar slant without antimicrobials

Acceptable: Sabouraud's dextrose agar slant or similar media without antimicrobials (eg, 7H11 agar slant, LJ, MGIT [7H9]

broth media)

Specimen Volume: Isolate

Collection Instructions: Organism must be in pure culture, actively growing.

Forms:

If not ordering electronically, complete, print, and send a Microbiology Test Request (T244) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)		
	Refrigerated		

Ask at Order Entry (AOE) Questions:

Test	Questi	Description	Туре	Report
ID	on ID			able
MM	Q00M0	Specimen Source (Required) and Organism Identification (Required unless	Plain	Yes
LNS	037	concurrent identification test is ordered)	Text	

Result Codes:

Result ID	esult ID Reporting Name		Unit	LOINC®
MMLNS	Susc, Aerobic Actinomycetes	Alphanumeric	mL	29577-4

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:



Test Definition: MMLNS

Antimicrobial Susceptibility, Nocardia species and other Aerobic Actinomycetes, Varies

N	\sim

CPT Code Information:

87186

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
MIC	Susceptibility, MIC	1	87186	No	No, (Bill Only)

Reference Values:

Interpretive criteria and reporting guidelines are followed using the Clinical Laboratory Standards Institute (CLSI) M24S document.