MAYO CLINIC LABORATORIES

Test Definition: FHSTW

Histamine, Whole Blood

Overview

Method Name
Quantitative Enzyme-Linked Immunosorbent Assay

NY State Available
Yes

Specimen

Specimen Type
WB Heparin

Shipping Instructions
CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.

Specimen Required
Specimen Type: Whole Blood
Collection Container/Tube:
Preferred: Green top (sodium heparin)
Acceptable: Green top (lithium heparin)
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL
Collection Instructions:
1. Draw blood in a green-top tube (sodium or lithium heparin).
2. Aliquot 1 mL well mixed whole blood into a plastic vial.
3. Freeze specimen.

Specimen Minimum Volume
0.5 mL

Reject Due To

Ambient or refrigerated samples	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
WB Heparin	Frozen	180 days	

Test Definition: FHSTW

Histamine, Whole Blood

Clinical & Interpretive

Reference Values
180 - 1800 nmol/L

Performance

PDF Report
No

Day(s) Performed
Monday, Thursday

Report Available
3 to 11 days

Performing Laboratory Location
ARUP Laboratories

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

CPT Code Information
83088

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FHSTW	Histamine, Whole Blood	46436-2

Result ID	Test Result Name	Result LOINC® Value
Z2753	Histamine, Whole Blood	46436-2