MAYO CLINIC LABORATORIES

Test Definition: FIMMC

Immunoglobulins, CSF Quantitative

Overview

Method Name
Quantitative Immunoturbidimetry

NY State Available
Yes

Specimen

Specimen Type
CSF

Specimen Required
Collection Container/Tube: Leak-proof container
Submission Container/Tube: Sterile vial
Specimen Volume: 1 mL
Collection Instructions:
1. Centrifuge spinal fluid to remove cellular material.
2. Transfer 1 mL into a plastic vial.
3. Send refrigerated.

Specimen Minimum Volume
0.6 mL

Reject Due To

Hemolysis:	Reject
Thawing:	Cold OK; Warm reject
Grossly bloody	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
CSF	Refrigerated (preferred)	30 days	
	Frozen	180 days	

Clinical & Interpretive

Clinical Information

Refer to https://ltd.aruplab.com/

Reference Values

Immunoglobulin M CSF: 0.0-0.7 mg/dL
Immunoglobulin G CSF: 0.0-6.0) mg/dL
Immunoglobulin A CSF: 0.0-0.7 mg/dL

Performance

PDF Report

No

Day(s) Performed

Wednesday, Saturday

Report Available

1 to 13 days

Performing Laboratory Location

ARUP Laboratories

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information

82784 x3

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FIMMC	Immunoglobulins, CSF Quantitative	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z2771	Immunoglobulin M CSF	2471-1
Z2772	Immunoglobulin G CSF	2464-6
Z2773	Immunoglobulin A CSF	2457-0