MAYO CLINIC LABORATORIES

Test Definition: FHADE

Haddock (Melanogrammus aeglefinus) IgE

Overview

Method Name
ImmunoCAP FEIA

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube:
Preferred: Red top
Acceptable: Serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 0.5 mL Serum
Collection Instructions:
1. Centrifuge and aliquot 0.5 mL of serum into a plastic vial.
2. Send refrigerate.

Specimen Minimum Volume
See Specimen Required

Reject Due To

Lipemia	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

Clinical & Interpretive

Clinical Information

Refer to www.eurofins-viracor.com/test-menu

Reference Values

<0.35 kU/L

Interpretation

Class	IgE (kU/L	Comment
0	<0.10	Negative
0/1	0.10 â€“ 0.34	Equivocal/Borderline
1	0.35 â€“ 0.69	Low Positive
2	0.70 â€“ 3.49	Moderate Positive
3	3.50 â€“ 17.49	High Positive
4	17.50 â€“ 49.99	Very High Positive
5	50.00 â€“ 99.99	Very High Positive
6	>99.99	Very High Positive

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

3 to 6 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

Test Definition: FHADE

Haddock (Melanogrammus aeglefinus) IgE

CPT Code Information

86003

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FHADE	Haddock IgE	7384-1

Result ID	Test Result Name	Result LOINC® Value
Z3112	Haddock IgE	7384-1
Z3113	Class	102358-9