MAYO CLINIC LABORATORIES

Test Definition: FADDS

Antidepressant Drug Screen, Qualitative

Overview

Method Name
Liquid Chromatography/Tandem Mass Spectrometry (LC/MS-MS)

NY State Available
Yes

Specimen

Specimen Type
Varies

Specimen Required
Submit only 1 of the following specimens:

Specimen Type: Plasma
Collection Container/Tube: Gray top (NaF/oxalate), green top (heparin), or lavender top (EDTA) (plasma gel/PST are **not acceptable**)
Submission Container/Tube: Plastic vial
Specimen Volume: 3 mL Plasma
Collection Instructions:
1. Centrifuge and aliquot 3 mL plasma into a plastic vial.
2. Send refrigerate.

Specimen Type: Serum
Collection Container/Tube: Red top (serum gel/SST are **not acceptable)**
Submission Container/Tube: Plastic vial
Specimen Volume: 3 mL Serum
Collection Instructions:
1. Centrifuge and aliquot 3 mL serum into a plastic vial.
2. Send refrigerate.

Specimen Type: Urine
Submission Container/Tube: Plastic preservative-free container
Specimen Volume: 3 mL
Collection Instructions:
1. Collect 3 mL from a random urine collection in a preservative-free plastic urine container.
2. Send refrigerated.

Specimen Minimum Volume

0.6 mL

Reject Due To

All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	
	Ambient	72 hours	
	Frozen	180 days	

Clinical & Interpretive

Clinical Information

Refer to https://medtox.labcorp.com/test-menu.

Interpretation

Antidepressant screen includes the analysis for:
Amitriptyline, clomipramine and desmethylclomipramine, cyclobenzaprine, desipramine, doxepin and desmethyldoxepin, fluoxetine and norfluoxetine, imipramine, maprotiline, nortriptyline, paroxetine, protriptyline, sertraline and desmethylsertraline, trimipramine

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

5 to 9 days

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information

80333
80337
80369

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FADDS	Antidepressant Drug Screen	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z3422	Specimen Type	31208-2
Z3423	Antidepressant Screen	Not Provided
Z3424	Amitriptyline	Not Provided
Z3425	Nortriptyline	Not Provided
Z3426	Imipramine	Not Provided
Z3427	Desipramine	Not Provided
Z3428	Sertraline	Not Provided
Z3429	Desmethylsertraline	Not Provided
Z3430	Fluoxetine	Not Provided
Z3431	Norfluoxetine	Not Provided
Z3432	Clomipramine	Not Provided
Z3433	Desmethylclomipramine	Not Provided
Z3434	Doxepin	Not Provided
Z3435	Desmethyldoxepin	Not Provided
Z3436	Maprotiline	Not Provided
Z3438	Trimipramine	Not Provided
Z3439	Paroxetine	Not Provided
Z3440	Cyclobenzaprine	Not Provided
Z3437	Protriptyline	Not Provided