MAYO CLINIC LABORATORIES

Test Definition: FWFHG

Venom W-F Hornet IgG

Overview

Method Name
Phadia ImmunoCAP

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube:
Preferred: Red top
Acceptable: Serum gel
Submission Container/tube: Plastic vial
Specimen Volume: 0.5 mL Serum
Collection Instructions:
1. Centrifuge and aliquot 0.5 mL of serum into a plastic vial.
2. Send refrigerated.

Specimen Minimum Volume
See Specimen Required

Reject Due To

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	7 days	
	Frozen	365 days	

Clinical & Interpretive

Clinical Information
Refer to www.eurofins-viracor.com/test-menu/

Reference Values

>3.5 mcg/mL

Venom IgG Reference Values: 1.0-3.5 mcg/mL
Low venom IgG indicating a significant reaction risk. 3.5-6.0 mcg/mL
Moderate level of venom IgG that may be associated from serious sting reactions(J Clin Immun.1983;6:172) >6.0 mcg/mL
Elevated venom IgG usually associated with protection from serious sting reactions(J Clin Immun.1983;6:172)
The interpretative guidelines have been adapted from Golden et al. (JACI. 1992;90:386-393).
The test method used is the Phadia ImmunoCAP IgG assay, which has been recalibrated to correlate with the venom IgG assay referenced above.

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

2 to 7 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

86001

LOINC® Information

MAYO CLINIC LABORATORIES

Test Definition: FWFHG

Venom W-F Hornet IgG

Test ID	Test Order Name	Order LOINC® Value
FWFHG	WF Hornet IgG	49745-3

Result ID	Test Result Name	Result LOINC® Value
Z3564	Venom W-F Hornet (D-maculata) IgG	49745-3