MAYO CLINIC LABORATORIES

Test Definition: FIL2M

Interleukin 2, Serum

Overview

Method Name
Quantitative Multiplex Bead Assay

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube: Serum gel or red top
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL Serum
Collection Instructions:
1. Within 2 hours of collection, centrifuge and aliquot 1 mL of serum into a plastic vial.
2. Freeze within 30 minutes and send frozen.
Note: **Critical frozen. Separate specimens must be submitted** when multiple tests are ordered.

Specimen Minimum Volume
Serum: 0.4 mL

Reject Due To

Contaminated specimens	Reject
Heat-inactivated specimens	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Frozen	365 days	

Clinical & Interpretive

Clinical Information

Refer to https://ltd.aruplab.com/

Reference Values
< or =2.1 pg/mL

Interpretation
Results are used to understand the pathophysiology of immune, infectious, or inflammatory disorders, or may be used for research purposes.

Cautions
Cytokine levels may demonstrate diurnal variation. For longitudinal comparison, it is recommended that cytokine levels be determined at the same time of day.

Performance

PDF Report
No

Day(s) Performed
Monday through Sunday

Report Available
3 to 8 days

Performing Laboratory Location
ARUP Laboratories

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by ARUP Laboratories.

CPT Code Information
83520

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value

FIL2M	Interleukin 2	33939-0

Result ID	Test Result Name	Result LOINC® Value
FIL2M	Interleukin 2	33939-0