MAYO CLINIC LABORATORIES

Test Definition: FDLS

Donath Landsteiner

Overview

Method Name
Hemolysis

NY State Available
Yes

Specimen

Specimen Type
Serum Red

Specimen Required
Collection Container/Tube: Red top (serum gel/SST are **not acceptable**)
Submission Container/Tube: Plastic vial
Specimen Volume: 3 mL serum
Collection Instructions:
1. Maintain specimen at 37 degrees C until after centrifugation.
2. Centrifuge and aliquot 3 mL of serum into a plastic vial.
3. Send refrigerated.

Specimen Minimum Volume
Serum: 2 mL

Reject Due To

Gross hemolysis	Reject
Gross lipemia	Reject
Gross icterus	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Ambient	72 hours	

Clinical & Interpretive

Clinical Information
Refer to https://ltd.aruplab.com/

Reference Values
Negative

Performance

PDF Report
No

Day(s) Performed
Monday through Friday

Report Available
1 to 7 days

Performing Laboratory Location
ARUP Laboratories

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information
86940
86941

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FDLS	Donath Landsteiner	18287-3

Result ID	Test Result Name	Result LOINC® Value
FDLS	Donath Landsteiner	18287-3