MAYO CLINIC LABORATORIES

Test Definition: FDGTX

Digitoxin, Serum

Overview

Method Name
Immunoassay (IA)

NY State Available
Yes

Specimen

Specimen Type
Serum Red

Specimen Required
*****Must submit one specimen per order. Specimens cannot be shared between multiple orders.*****

Collection Container/Tube: Red top (serum gel/SST are **not acceptable**)
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL Serum
Collection Instructions:
1. Promptly centrifuge and aliquot serum into a plastic vial.
2. Send refrigerate.

Specimen Minimum Volume
Serum: 0.4 mL

Reject Due To
All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Frozen	90 days	

Clinical & Interpretive

Clinical Information
Refer to www.nmslabs.com/test-catalog.

Reference Values

Reporting limit determined each analysis.

None detected
Therapeutic range: 10-30 ng/mL

Performance

PDF Report

Referral

Day(s) Performed

Monday through Sunday

Report Available

9 to 11 days

Performing Laboratory Location

NMS Labs

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information

80299

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FDGTX	Digitoxin	3559-2

Result ID	Test Result Name	Result LOINC® Value
FDGTX	Digitoxin	3559-2