MAYO CLINIC LABORATORIES

Test Definition: FACES

Acetoacetate, Serum

Overview

Method Name
Gas Chromatography (GC)

NY State Available
Yes

Specimen

Specimen Type
Serum Red

Specimen Required
*****Must submit one specimen per order. Specimens cannot be shared between multiple orders.*****

Collection Container/Tube: Red top (Serum gel/SST are **not acceptable**)
Submission Container/Tube: Plastic vial
Specimen Volume: 3 mL Serum
Collections Instructions:
1. Centrifuge and aliquot 3 mL of serum into a plastic vial.
2. Send frozen.

Specimen Minimum Volume
Serum: 1.2 mL

Reject Due To
All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Frozen (preferred)	4 days	
	Refrigerated	4 days	

Clinical & Interpretive

Clinical Information
Refer to www.nmslabs.com/test-catalog

Reference Values

Reporting limit determined each analysis.

Acetoacetate: None Detected
Normal range for adults: 5-30 mcg/mL

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

5 to 11 days

Specimen Retention Time

2 weeks

Performing Laboratory Location

NMS Labs

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

82010

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FACES	Acetoacetate, Serum	Not Provided

Test Definition: FACES

Acetoacetate, Serum

Result ID	Test Result Name	Result LOINC® Value
Z3524	Acetoacetate	Not Provided